Description

Tetric EvoFlow is a flowable, light-curing, radiopaque nano-hybrid composite for the restorative therapy and the cementation of ceramic and composite restorations. Tetric EvoFlow cures with light in the wavelength range of 400–500 nm (blue light).

Shades

Tetric EvoFlow is available in several shades.

Composition

The monomer matrix is composed of dimethacrylates (38 % wt). The fillers comprise barium glass, ytterbium trifluoride, highly dispersed silicon dioxide, mixed oxide and copolymer (62 % wt). Additives, catalysts, stabilizers and pigments are additional contents (< 1 % wt).

The total content of inorganic fillers is 57.5 % wt or 30.7 % vol. The particle sizes of the inorganic fillers range between 40 nm (0.04 µm) and 3000 nm (3 µm). The mean particle size is 550 nm (0.55 µm).

Indication

- As initial layer under Class I to V restorations
- Anterior restorations (Classes III, IV)
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Small restorations of all types
- Extended fissure sealing
- Splinting of mobile teeth
- Blocking out of undercuts
- Adhesive cementation of light-transmissive indirect composite and ceramic restorations
- Repair of composite and ceramic veneers

Contraindication

The use of Tetric EvoFlow is contraindicated:

- if an adequately dry working field cannot be established or if the stipulated technique cannot be applied;
- if the patient is known to be allergic to any of the ingredients in Tetric EvoFlow.

Side effects

In individual cases, components of Tetric EvoFlow may lead to sensitization. Tetric EvoFlow should not be used in such cases.

To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (apply a calcium hydroxide-based preparation to areas in close proximity to the pulp).

Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of composite materials. Consequently, application of such materials together with Tetric EvoFlow must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Tetric EvoFlow®

Instructions for Use	Bruksanvisning	– Light-curing, flowable resin-based dental restorative material	– Material de restauración dental en base a resina, fotopolimerizable y fluido
Gebrauchsinformation	Bruksanvisning	– Lichthärtendes, fließfähiges, zahnärztliches Füllungscomposit	– Composito de restauração dental fluido e fotopolimerizável
Mode d'emploi	Käyttöohjeet	– Composite de restauration fluide, photopolymérisable	– Ljushärdande flytande resinbaserat dentalt fyllnadsmaterial
Istruzioni d'uso	Bruksanvisning	– Composito fotopolimerizzabile fluido ad uso odontoiatrico	
Instrucciones de uso	Productinformatie		
Instruções de Uso	Οδηγίες Χρήσεως		
	Kullanma Talimatı		

Tetric EvoFlow®

Navodilo za uporabo	Informazione per �rdorim	– Teko�i smolasti dentalni obnovevitveni material, ki se strjuje s svetlobo	– Sv�tlem vytvrzovan�, tekut�, dent�ln� material pro restorativn� l�ebu zalo�en� na pryskyřici
Upute za uporabu	Instru�ioni de utilizare	– Proto�ni stomatološki kompozit za restauraciju sa svjetlosnim stvrdnjavanjem	– Tekut�, svetlom vytvrzovan� material na b�ze �ivice pre dent�lnu rekonstruk�n� lie�bu.
N�vod k pou�it�	Инструкция по использованию		
N�vod na pou�itie	Информ�ция щодо використання		
Haszn�lati utasit�s	Касутамисјуhend		
Instrukcja stosowania	Lietošanas inform�cija		
Упутство за употребу	Vartojimo instrukcija		
Упатство за употреба			
Инструкции за употреба			

Application in restorative therapy

1. Select the shade (see shade guide).
2. Provide appropriate isolation, best with a rubber dam (e.g. OptraDam® Plus).
3. Prepare the cavity according to the requirements of the adhesive technique.
4. Clean the cavity with water spray.
5. Air dry the cavity.
6. Do not apply a base material when using an enamel/dentin bonding agent. Selectively cover very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently use a pressure-resistant cement (e.g. a glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.
7. Apply matrix (e.g. OptraMatrix®) / interdental wedge
8. Condition and apply the bonding agent according to the instructions for use of the product in use. We recommend using Syntac® (with phosphoric acid etching) or ExciTE® F (with phosphoric acid etching) or the AdheSE® One F self-etching adhesive.
9. Apply Tetric EvoFlow in layers of max. 2 mm or 1.5 mm (Dentin shades, Bleach XL). Polymerize each layer according to the table shown below. Hold the light emission window as closely as possible to the surface of the restorative material. When using a metal matrix, additionally polymerize from the buccal or lingual/palatal after removing the matrix.
10. Tetric EvoFlow curing times: see Table 1 (at the end of the Instructions for use).
11. Finish the restoration with suitable finishers (e.g. Astropol® F or the single-step OptraPol® Next Generation) or fine diamonds.
12. Check the occlusion.
13. Polish with silicone polishers (e.g. Astropol P/HP, Astrobrush® or the single-step OptraPol® Next Generation) as well as with polishing discs and polishing strips.

Application for the cementation of ceramic or composite restorations

Tetric EvoFlow is to be used exclusively in light-curing procedures. Consequently, the composite is only suitable for the cementation of restorations that are sufficiently translucent to permit enough polymerization light to pass through the restoration and reach Tetric EvoFlow in order that the composite can cure completely. Before being cemented in place, the restoration has to be pre-treated according to the instructions of the manufacturer.

Prepare the cavity as follows:

1. Appropriate isolation, i.e. with a rubber dam or OptraDam Plus, is required.
2. Remove the temporary restoration and any residue of the temporary luting material and clean the cavity with a prophylactic paste (e.g. Proxyl®); subsequently, rinse with water spray and dry the cavity.
3. Condition and apply the bonding agent according to the instructions for use of the product in use. If necessary, protect the adjacent tooth with a matrix band during conditioning.
4. Apply Tetric EvoFlow directly to the restoration; insert the restoration, and remove excess material. To prevent the formation of an oxygen inhibited layer, an oxygen-tight glycerine gel (e.g. Liquid Strip) may be applied along the restoration margins. Make sure that the glycerine gel is not mixed with Tetric EvoFlow!
5. Next, polymerize Tetric EvoFlow in segments starting with the proximal margin:
 - If a curing unit with a light intensity of min. 800 mW/cm² is used (e.g. bluephase®), light-cure for 10 s per mm of ceramic and segment.
 - If a curing unit with a light intensity of 500 to 800 mW/cm² is used, light-cure for 20 s per mm of ceramic and segment.Remove the gel with water spray after polymerization.
6. After polymerization, remove excess material with suitable finishers (e.g. Astropol F or single-step OptraPol Next Generation) or fine diamonds. Remove proximal excess with diamond or tungsten carbide finishers or finishing strips. Check the occlusion and articulation and adjust if necessary. Use silicone polishers (e.g. Astropol P/HP, OptraPol Next Generation) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Tetric EvoFlow should have room temperature when applied.
2. If necessary, fresh Tetric EvoFlow can be directly applied to polymerized material. If the material has already been polished, it must first be roughened before fresh Tetric EvoFlow is applied.

Warning

- Unpolymerized Tetric EvoFlow should not come in contact with skin/mucous membrane or eyes. Unpolymerized Tetric EvoFlow may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Storage temperature: 2–28 °C / 36–82 °F.
- Close syringes/Cavifils immediately after usage. Exposure to light causes premature polymerization.
- Shelf life: see note on Cavifils, syringes and packages.
- Do not use Tetric EvoFlow after the indicated date of expiration.
- Do not use oxidizing disinfectants to disinfect syringes and Cavifils.
-  If Tetric EvoFlow is applied from the Cavifil directly into the oral cavity of the patient, we recommend using the Cavifil only once due to hygienic reasons (prevention of cross-contamination between patients). The same applies to the application tips of syringes.

Keep the material out of the reach of children.

For use in dentistry only.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Table 1

Program \ Unit	bluephase C8	bluephase	bluephase 20i	bluephase style
Turbo	–	–	5 s	–
High Power	15 s	10 s	10 s	10 s
Soft Start	20 s	15 s	15 s	–

Light Intensity	Exposure Time
$\geq 500 \text{ mW/cm}^2$	20 s
$\geq 1000 \text{ mW/cm}^2$	10 s

Ivoclar Vivadent AG

Bendererstrasse 2, FL-9494 Schaan, Liechtenstein
Tel. +423 235 35 35, Fax +423 235 33 60,
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia
Tel. +61 3 979 595 99, Fax +61 3 979 596 45
www.ivoclarvivadent.com.au

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723, Centro Empresarial Tamboré,
CEP 06460-110 Barueri – SP, Brazil
Tel. +55 11 2424 7400, Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

2785 Skymark Avenue, Unit 1, Mississauga, Ontario L4W 4Y3,
Canada, Tel. +1 905 238 5700, Fax +1 905 238 5711
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Rm 603 Kuen Yang, International Business Plaza, No. 798 Zhao Jia Bang Road,
Shanghai 200030, China
Tel. +86 21 5456 0776, Fax +86 21 6445 1561
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia
Tel. +57 1 627 33 99, Fax +57 1 633 16 63
www.ivoclarvivadent.com

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France
Tel. +33 450 88 64 00, Fax +33 450 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany
Tel. +49 (0) 79 61 / 8 89-0, Fax +49 (0) 79 61 / 63 26
www.ivoclarvivadent.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate, Veera Desai Road, Andheri (West),
Mumbai, 400 053, India
Tel. +91 (22) 2673 0302, Fax +91 (22) 2673 0301

Ivoclar Vivadent s.r.l.

Via Isonzo 67/69, 40033 Casalecchio di Reno (BO), Italy
Tel. +39 051 611 35 55, Fax +39 051 611 35 65
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.

1-28-24-4F Hongo, Bunkyo-ku, Tokyo 113-0033, Japan
Tel. +81 3 6903 3535, Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37, Seocho-dong, Seocho-gu, Seoul 137-855,
Republic of Korea, Tel. +82 (2) 536 0714, Fax +82 (2) 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863, Piso 14, Col. Napoles, 03810 México, D.F., Mexico
Tel. +52 (55) 5062-1000, Fax +52 (55) 5062-1029
www.ivoclarvivadent.com.mx

Ivoclar Vivadent Ltd.

12 Omega St, Albany, PO Box 5243 Wellesley St, Auckland, New Zealand
Tel. +64 9 914 9999, Fax +64 9 814 9990
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78, PL-00175 Warszawa, Poland
Tel. +48 22 635 54 96, Fax +48 22 635 54 69
www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Derbenevskaja Nabereshnaya 11, Geb. W, 115114 Moscow, Russia
Tel. +7 495 913 66 19, Fax +7 495 913 66 15
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204,
P.O. Box 300146, Riyadh 11372, Saudi Arabia
Tel. +966 1 293 83 45, Fax +966 1 293 83 44
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

171 Chin Swee Road, #02-01 San Centre, Singapore 169877
Tel. +65 6535 6775, Fax +65 6535 4991
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

c/ Emilio Muñoz Nº 15, Entrada c/ Albarracin, E-28037 Madrid, Spain
Tel. + 34 91 375 78 20, Fax + 34 91 375 78 38
www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, S-169 56 Solna, Sweden
Tel. +46 (0) 8 514 93 930, Fax +46 (0) 8 514 93 940
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

Ahi Evran Caddesi No 1, Polaris Is Merkezi Kat: 7, 80670 Maslak, Istanbul, Turkey, Tel.
+90 212 346 04 04, Fax +90 212 346 04 24
www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building, Feldspar Close, Warrens Business Park,
Enderby, Leicester LE19 4SE, United Kingdom
Tel. +44 116 284 78 80, Fax +44 116 284 78 81
www.ivoclarvivadent.com

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA
Tel. +1 800 533 6825, Fax +1 716 691 2285
www.ivoclarvivadent.com