

Coe® Tray Cleaner

Dustless Formulation

For Removing Alginate, ZOE, Waxes and Modeling Compound

Directions:

1. Stir two capfuls (2 oz.) of powder into one quart of warm water until completely dissolved.
2. Strip bulk of material from instruments or trays and completely immerse in solution.
3. Immerse the tray in a solution of tray cleaner and warm water for 15-30 minutes. Hot solution removes residue in a few minutes.
4. Stubborn particles may be removed with slight brushing. Rinse instruments and tray with running water before processing or storing.

Caution: Clean nickel and chrome instruments and tray in a separate solution from that used for aluminum instruments. Do not leave instruments in solution longer than is necessary for cleaning. Prolonged immersion may cause some corrosion, particularly on worn areas.

Contains: Sodium Carbonate

- Causes burns
- Irritating to eyes and respiratory system
- Do not breath dust
- Avoid contact with skin and eyes
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
- Wear suitable protective clothing and gloves



GC America INC.
Alsip, IL 60803

Made in U.S.A.