ENGLISH
<b>Temp-Bond™</b> - Zinc Oxide/Eugenol Cement
<b>Temp-Bond™ NE</b> - Zinc Oxide/Non-Eugenol Cement
<b>Temp-Bond™ Clear with Triclosan</b> - Dual Curable Transparent Resin Cement with Triclosan
<b>Temporary Cements</b>
<b>DESCRIPTION</b>
Temp-Bond temporary cement product series consists of Temp-Bond, Temp-Bond NE and Temp-Bond Clear with Triclosan. Each is designed to be suitable for various applications such as cementing temporary crowns, bridges, inlays, onlays and splints. It has excellent flow to permit the restoration to be easily and completely seated. It is strong enough to withstand the stresses of mastication, yet permits easy removal of the restoration when this is desired. Temp-Bond and Temp-Bond NE are available in traditional tubes, Unidose™ single-use foil pouches, and automix dual barrel syringes. Temp-Bond Clear with Triclosan is available in an automix dual barrel syringe.

**DIRECTIONS FOR USE**

*Work Time and Set Time:*

At the ideal mixing conditions of ambient temperature, the following work time and setting time from start of mixing are obtained:

	Work Time (mm:ss)	Set Time (mm:ss)
Temp-Bond	≥ 1:30	≤ 7:00
Temp-Bond NE	≥ 1:30	≤ 7:00
Temp-Bond Clear with Triclosan	≥ 1:30	≤ 5:30

These times may vary based on storage conditions, temperature, humidity, and age of product.

*Mixing:*

Tube: Extrude equal lengths of Base and Accelerator onto the mixing pad provided.

The length to be extruded will depend on the size and type of restoration to be cemented. Replace caps tightly on tubes after use. Thoroughly mix the pastes for approximately 30 seconds.

**Temporary Cements**

**DESCRIPTION**

Temp-Bond temporary cement product series consists of Temp-Bond, Temp-Bond NE and Temp-Bond Clear with Triclosan. Each is designed to be suitable for various applications such as cementing temporary crowns, bridges, inlays, onlays and splints. It has excellent flow to permit the restoration to be easily and completely seated. It is strong enough to withstand the stresses of mastication, yet permits easy removal of the restoration when this is desired. Temp-Bond and Temp-Bond NE are available in traditional tubes, Unidose™ single-use foil pouches, and automix dual barrel syringes. Temp-Bond Clear with Triclosan is available in an automix dual barrel syringe.

**DIRECTIONS FOR USE**

*Work Time and Set Time:*

At the ideal mixing conditions of ambient temperature, the following work time and setting time from start of mixing are obtained:

	Work Time (mm:ss)	Set Time (mm:ss)
Temp-Bond	≥ 1:30	≤ 7:00
Temp-Bond NE	≥ 1:30	≤ 7:00
Temp-Bond Clear with Triclosan	≥ 1:30	≤ 5:30

These times may vary based on storage conditions, temperature, humidity, and age of product.

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