

# DURELON®/DURELON FAST SET TECHNIQUE



## STEP 1:

Remove provisional restoration.



## STEP 2:

Thoroughly clean preparation with water.




## STEP 3:

Gently dry preparation. Do not over dry.




## STEP 4:

Dry the restoration with air.




## STEP 5:

Dispense one spoonful of powder onto mixing pad.




## STEP 6:

Dispense liquid.




## STEP 7:

Mix powder and liquid.




## STEP 8:

Place a thin, even coat of cement on inner surface of restoration. Do not use the mixture if it has lost its gloss.



## STEP 9:

Seat the restoration without pressure.



## STEP 10:

Allow cement to set, then remove excess.

# DURELON®/DURELON FAST SET

**C**arboxylate cement in the traditional powder and liquid form for cementing inlays, crowns, bridges, and more.

## INDICATIONS

- Inlay cementation
- Crown and bridge cementation
- Cement bases

## ADVANTAGES

- Low film thickness ensures optimized fit and marginal integrity
- Fluoride-releasing
- Chemically adheres to both enamel and dentin for reduced microleakage
- Biocompatible for sensitive patients
- Relatively neutral pH value for increased patient acceptance
- Less soluble than zinc phosphate or ZOE cements
- Also available in a fast set version that sets twice as fast with the same physical properties



## PACKAGING

### Type of Package:

### Item Number

#### DURELON Polycarboxylate Cement

INTRO PACK - Powder (20g), Liquid (14g)	.038019
TRIPLE PACK - Powder (60g), Liquid (40g)	.038233
TRIPLE LIQUID - (40g)	.038218
TRIPLE POWDER - Regular (60g)	.038206
TRIPLE POWDER - Fast Set (60g)	.038507

## MIXING RATIOS/WORKING & SETTING TIMES

Application	Powder/Liquid Ratio	Working Time (including mixing time)	Setting Time (from start of mix)
Cementation	1 scoop powder to 3 drops liquid	2½ minutes	10 minutes
Relining	1 scoop powder to 2 drops liquid	2 minutes	8 minutes

*Note: Fast set powder sets in half the time as regular set powder.*

TO PLACE AN ORDER, CALL YOUR  
AUTHORIZED ESPE DISTRIBUTOR  
OR CALL ESPE AT 1-800-344-8235.