Venus[®] flow

- GE** Gebrauchsanweisung
- GB** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- NL** Gebruiksaanwijzing

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Manufacturer:
Kulzer GmbH
Leipzigstraße 2
63450 Hanau (Germany)

Made in Germany

Distributed in USA / Canada exclusively by:
Kulzer, LLC
4315 South Lafayette Blvd.
South Bend, IN 46614-2517
1-800-431-1785

Federal law restricts this device to sale by or on the order of a dentist.

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Venus[®] flow

Venus Flow is a flowable, light curing, radiopaque hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations. Due to its low viscosity Venus Flow provides excellent coverage of the tooth structure.

Composition:
Venus Flow is based on multifunctional methacrylate monomers (EBADMA / TEGDMA) and contains approximately 62 % by weight or 38 % by volume inorganic fillers such as Ba-Al-F silicate glass and SiO₂. The filler particle sizes is between 0.005 µm and 5 µm.

Shades: A1, A2, A3, A3.5, A4, B2, B3, OA2, SBI, SB2, SBO, T2, HKA2.5, Baseline

Applications:

- Fissure sealing
- Enlarged fissure sealing
- Class V fillings
- Minimally invasive Class I and II fillings in areas not subjected to masticatory forces
- Minimally invasive Class III fillings
- Smaller contour and shade adjustment on the enamel and dentine
- Small surface repairs to direct and indirect restorations in combination with a suitable bonding agent
- Splinting of teeth
- Cavity lining
- Brackets retention

Kind of application:
Prepare the tooth structure for treatment by cleaning it with a fluoride free polishing paste. Before treatment, while the teeth are still moist, select the Venus Flow shade using the Venus 2Layer shade guide.

1. Remove any caries and complete minimally invasive preparation of the tooth structure according to the general guidelines for adhesive fillings. Then clean and dry the cavity. The use of a rubber dam is recommended.

2. A quality adhesive system (e.g. Bond[®] Self Etch, GLUMA[®] Comfort[®] Bond + Desensitizer, GLUMA[®] Solid Bond) should be applied to the enamel and dentin in accordance with the manufacturer's instructions before inserting the filling material. Avoid any contamination of conditioned or adhesive-treated surfaces with moisture, blood or saliva.

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Instructions for use **GB**

3. Apply Venus Flow in thin layers (max. 2 mm, Baseline max. 1 mm) to the cavity. When applying Venus Flow ensure that it flows slowly and evenly to prevent air bubbles – apply a gentle stream of air to facilitate maximum marginal contact. If necessary, remove air bubbles with a probe before polymerisation.
4. Polymerise each layer according to the table below using a suitable light-curing unit (e.g. Translux[®]) for 20 or 40 sec. A dispersion layer forms during polymerisation. This layer should not be touched or removed if other composite layers are to be applied.

Shades	Time
A1, A2, A3, A3.5, B2, B3, SBI, SB2, T2, HKA2.5	20 s
A4, OA2, SBO, Baseline	40 s

5. Venus Flow can be prepared and polished immediately following polymerisation. Finishing diamonds, tungsten carbide instruments, flexible discs, silicone polishers and polishing brushes are suitable for preparing. During long restorative procedures, the surgery light should be temporarily directed away from the treatment area to prevent premature initiation of the polymerisation reaction.

Note:

- A Kulzer Translux[®] light unit or a light unit with comparable intensity (min. 400–500 mW/cm²) should be used. The same polymerisation time is required with an LED light unit. If a stronger polymerisation unit (e.g. a plasma light unit with a greater output than 1200 mW/cm²) is used, the polymerisation time can be reduced.
- The polymerisation times given in the instructions for use are the minimum times required to attain optimum polymerisation.
- During polymerisation the light guide should be placed as close as possible to the surface being polymerised.
- The tips supplied for hygienic application are intended for single-use and should be discarded after use. After removing the material from the syringe take the pressure off the plunger to prevent continued flow. Replace the cap of the syringe after use.
- With small surfaces in caries-free exposed pulp cavities, the exposed area should be directly capped before beginning the adhesive treatment (e.g. a combination of CaOH and glass ionomer cement). As small an area as possible of the dentine is wetted with the curing light (e.g. glass ionomer cement).
- Do not use Venus Flow with eugenol cements, as these could inhibit polymerisation.

Avoid contact with the skin or eyes. On contact with the eyes, rinse with plenty of water and consult a doctor if necessary.

Do not use Venus Flow after the expiry date or after 25°C/77°F. Avoid direct sunlight. This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer

For use by dental personnel only: Syringes 1 x 1,8 g

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