GC FujiCEM™ 2

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Cementation of metal-based inlays, onlays, crowns and bridges
2. Cementation of resin inlays, onlays, crowns and bridges
3. Cementation of all ceramic inlays
4. Cementation of high strength (e.g. zirconia based, lithium disilicate, etc.) all ceramic inlays, onlays, crowns and bridges
5. Cementation of metal, ceramic and fiber posts

CONTRAINDICATIONS

1. Pulp capping.
2. In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

IDENTIFICATION OF PARTS (Fig. 1)

1. Paste Pak Cartridge SL (Slide and Lock System)
1. Cover 2. Inner plunger 3. Cartridge tip 4. Sleeve
- B. Paste Pak Dispenser
5. Piston 6. Adjusting guide 7. Lever 8. Sliding block 9. Piston release
- C. GC FujiCEM 2 Mixing Tip SL
- GC FujiCEM 2 Mixing Tip SL with intraoral tip

DIRECTIONS FOR USE

Paste / Pasta ratio (g/g)	2.0 / 1.0
Mixing time (sec.)	10'-15"
Working time (min., sec.) (at 23°C) (from start of mixing)	2'15"
Final finishing commencing time (min., sec.) (after seating the restoration)	4'30"

PASTE PAK CARTRIDGE SL LOADING

1. Make sure that the piston is completely retracted into the dispenser (Fig. 2).
Note:
To retract the piston, push the piston release and at the same time pull the sliding block (on the underside of the dispenser) with fingers.
2. Load the cartridge into the Paste Pak Dispenser (Fig. 3). Ensure the ▲ mark on the cartridge lines up with the ▲ mark on the Paste Pak Dispenser as illustrated. Then turn the cartridge fully to the right position.
3. After cartridge loading, push the piston release forward until it stops (Fig. 4).
a) Push the piston release forward.
b) Do not push the upper slide.
Note:
To remove the cartridge, retract the piston fully and turn the cartridge back to the marked positions.

A. CEMENTATION TECHNIQUE – AUTOMIXING WITH GC FujiCEM 2 MIXING TIP SL

1. TOOTH PREPARATION
a) Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
b) Clean the prepared tooth with pumice and water.
Note:
1) NaClO, EDTA or Chlorhexidine can be used to clean the prepared tooth. Do not use substances like H₂O₂ (hydrogen peroxide) and other disinfectants as these may lower bond strength.
2) Use GC Fuji PLUS Conditioner on the prepared tooth surface for 20 seconds for optimal adhesion.
3) When increased retention is required for ceramic prostheses use a ceramic primer such as CERAMIC PRIMER II (GC) and follow manufacturer's instructions.
- 2) Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. RESTORATION PREPARATION
Make sure that the restoration is pretreated and handled according to the manufacturers instructions.

3. DISPENSING AND AUTOMIXING
a) Remove the cartridge cover.
b) Before every use, check the two openings of the cartridge tip and confirm the pastes are at the same level to ensure even flow from the cartridge.
Note:
Before dispensing the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge.

- a) Attach the GC FujiCEM 2 Mixing Tip SL on the cartridge. Push the mixing tip completely until the end (Fig. 5-1). Hold the sides of the sleeve and slide it forward to lock the mixing tip (Fig. 5-2).
Note:
1) Make sure GC FujiCEM 2 Mixing Tip SL is securely attached to the cartridge. If not, remove the GC FujiCEM 2 Mixing Tip SL and reattach.
2) If the material has been refrigerated, leave the cartridge in room temperature for 10 minutes.

- d) Set the adjusting guide to the front of the lever (Fig. 6).
Note:
Before dispensing the material into the appliances, bleed a small amount of the mixed pastes on a pad.
e) Slowly depress the lever to dispense the material into the appliance.

- f) After dispensing, release the lever, slide the sleeve backward and remove the GC FujiCEM 2 Mixing Tip SL (Fig. 7). Clean the cartridge tip and always replace the cartridge cover after use.
Note:
1) Do not leave the used GC FujiCEM 2 Mixing Tip SL on the cartridge for a long period of time.
2) The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser. As the cartridge cylinders are translucent, the residual amount can also be confirmed by the position of the inner plunger seen through the cylinders. When only a small amount of paste remains in the cartridge, the paste may not be dispensed from the GC FujiCEM 2 Mixing Tip SL.
3) To avoid damage, do not drop the dispenser with the cartridge and mixing tip.

4. CEMENTATION
a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 15 seconds from start of mixing at 23°C (74°F). Higher temperature will shorten the working time.
b) Maintain moderate pressure and make sure the restoration remains in place. Start removing excess cement with a scaler when the excess cement feels bubbly.
c) Finishing can be started 4 minutes 30 seconds after seating the restoration.

- B. CEMENTATION TECHNIQUE – HANDMIXING

1. TOOTH PREPARATION
Refer to A. 1.

2. RESTORATION PREPARATION
Refer to A. 2.

3. DISPENSING
a) Remove the cartridge cover.
b) Depress the lever to dispense the required amounts of paste onto the mixing pad.

- 1) Before dispensing the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge.
2) The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed (Fig. 8).
1. Position of adjusting guide for one inlay.
2. Position of adjusting guide for one crown.

- c) While depressing the lever (Fig. 9-1), level and cut off the extruded material from the cartridge tips by moving the cartridge and dispenser to an upright position on the mixing pad (Fig. 9-2).
- d) When the lever is released after dispensing, residual pastes at the cartridge tips are retracted into the cartridge. Clean the cartridge tip and replace the cartridge cover.
Note:
The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser. As the cartridge cylinders are translucent, the residual amount can also be confirmed by the position of the inner plunger seen through the cylinders.

4. MIXING
After dispensing, spread the two pastes thinly on the mixing pad using a plastic spatula. Mix thoroughly, with lapping strokes, for 10 seconds. Be careful not to incorporate air bubbles.

- Note:
1) Larger amount of paste requires a longer mixing time (15 seconds).
2) Do not mix GC FujiCEM 2 pastes with other materials.

5. CEMENTATION
Refer to A. 4. Refrigeration of the cartridge is recommended to extend the working time.

- CLEANING / STERILIZATION OF PASTE PAK DISPENSER
For cleaning and sterilization of PASTE PAK DISPENSER, refer to the respective instructions for use.

- SHADE
Light Yellow

- STORAGE
Recommended for optimal performance, store at room temperature away from heat, moisture and direct sunlight (4-25°C / 39-77°F).
Shelf life: 2 years from date of manufacture.

- CLEANING AND DISINFECTING RECOMMENDATION
MULTI-USE DELIVERY SYSTEMS. To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Do not use if device is damaged.

- DO NOT IMMERSIVE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

- PACKAGES
1. GC FujiCEM 2 Automatic
Paste Pak Cartridge SL (13.3g / 7.2mL) (2), GC FujiCEM 2 Mixing Tip SL (42)

2. GC FujiCEM 2 Refill
Paste Pak Cartridge SL (13.3g / 7.2mL) (2), with mixing pad (No. 22)

3. GC FujiCEM 2 Starter Kit
Paste Pak Cartridge SL (13.3g / 7.2 mL) (1), GC FujiCEM 2 Mixing Tip SL (20), GC FujiCEM 2 Dispenser (1)
4. Paste Pak Dispenser (1)

- CAUTION
Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. In rare cases metal may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of Paste Pak Dispenser (nickel-plated aluminum) and refer to a physician.

4. Store Paste Pak Dispenser away from chemicals such as iodine tincture containing halogen element. Otherwise the dispenser may discolor.

5. GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.

6. Do not use GC FujiCEM 2 in combination with eugenol containing materials as eugenol may hinder GC FujiCEM 2 from setting or bonding properly.
7. This product is not indicated for filling or core build-up.

8. The GC FujiCEM 2 Mixing Tip SL (including the intraoral tip) is for single use only. To prevent cross contamination between patients, do not reuse the GC FujiCEM 2 Mixing Tip SL. The GC FujiCEM 2 Mixing Tip SL cannot be sterilized in an autoclave or chemically.

- Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available for the Americas at: <http://www.gcamerica.com>

- They can also be obtained from your supplier.
Last Revised: 08/2017

- R.Only
MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.

- DISTRIBUTED BY
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
www.gcamerica.com

- GC SOUTH AMERICA
Rua Heliópolis, 399, Santana - São Paulo, SP - Brasil
CEP: 02022-031 - TEL: +55-11-2323-0865
CNFJ: 08.279.999.0001-61
RESP: TEC. Mayara De Santos Ribeiro - CRO/SP: 105.992
PRINTED IN THE U.S.A.

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