Acroseal

Root canal filling material by sealing gutta percha points

Composition

Base: Bismuth carbonate, esterified colophane, TCD Diamine-Aerosil.
Catalyst: Bismuth carbonate, calcium hydroxide, diglycidyl ether bisphenol A, yellow iron oxide.

Properties

ACROSEAL is a calcium hydroxide-based cement, made up of two components – a BASE paste and a CATALYST paste – and designed for the permanent filling of root canals by sealing up gutta percha points. ACROSEAL contains no eugenol.

Indications

Root canal filling by sealing gutta percha points.

Contraindications

Hypersensitivity to any component of the formula.

Instructions for use

Isolate the tooth to be filled.

- **Tubes:**
Lay down two equal lengths of BASE and CATALYST squeezed out from both tubes, on a mixing pad.
Spatulate the two pastes until you obtain homogeneous consistency.
- **“Automix” syringe:**
Attach a mixing tip to the double-syringe and lay down the required amount of paste onto a mixing pad.
Coat the gutta percha point(s) with ACROSEAL and condense the point(s).

Special warnings and precautions for use

As with any root canal filling material, only a perfectly shaped and disinfected canal can be filled. Determine the working length of the canals before filling, to prevent apical extrusion.

Special risks:

Base/ Catalyst: Corrosive. Causes burns. May cause sensitization by skin contact.

Safety advice:

Avoid contact with eyes and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water and soap.

Remarks

- When opening the tube of BASE, eliminate, if necessary, the supernatant liquid layer
- The setting time of ACROSEAL, in vivo, ranges from 3 and 12 hours and varies according to hygrometry
- Spatulas, mixing pads or any other instrument should be cleaned with alcohol or acetone, immediately after use

Storage

Do not store above 77°F (25°C).
Tightly close the tubes after each use.

Presentation

- | | |
|------------------------------------------------------|------------------------------------------------------------|
| Box containing:
• 1 x 0.3 oz (8.5 g) tube of BASE | Box containing:
• 2 x 0.3 oz (8.6 g) “Automix” syringes |
| • 1 x 0.34 oz (9.5 g) tube of CATALYST | • 15 mixing tips |
| • 1 cross-ruled mixing pad | • 15 intraoral tips |
| | • 1 mixing pad |

Rx only

For professional dental use only

Distributed in USA by: SEPTODONT
Louisville CO 80027 - USA
Distributed in Canada by: SEPTODONT
Cambridge, ON N1R 6X3 - CANADA
Made by SEPTODONT
94100 Saint-Maur-des-Fossés - France

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