**(en) Polyether Impression Material**

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Instructions for Use  
Mode d'emploi  
Información de uso

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Dental Products  
Carl-Schurz-Str. 1  
41453 Neuss – Germany

3M ESPE Customer Care/MSDS Information:  
U.S.A. 1-800-634-2249 and Canada 1-888-363-3685.

3M, ESPE, Garant, Impregum and Permadyne are trademarks of 3M or 3M Deutschland GmbH.  
Used under license in Canada. © 2016, 3M. All rights reserved.

44000018162/02

**(en) ENGLISH**

**Product Description**

Impregum™ F is a medium-bodied consistency polyether impression material for hand mixing. The mixing ratio is 7 volumes base paste : 1 volume catalyst.

For details on Polyether Adhesive, Elastomer syringe, Permadyne™ Garant™ 2:1 and Permadyne™ light-bodied consistency please refer to the corresponding instructions for Use. These instructions for use should not be discarded for the duration of product use.

**Indications**

- Impression of inlay, onlay, crown, bridge and veneer preparations
- Functional impressions
- Fixation impressions
- Implant impressions

**Precautionary Measures**

**For Patients**

The base paste of this product contains substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an adverse reaction is observed in a patient discontinue use of the product in this patient and seek medical attention as needed.

Known allergens are 1-Dodecylimidazole and Carthusian Essence (different flavors).

**For Dental Personnel**

The base paste of this product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended.

3M MSDSs can be obtained from [www.mmm.com](http://www.mmm.com) or contact your local subsidiary.

**Preparation**

Block-out undercuts or areas where gingival recession exists to prevent the impression material from “locking” onto the tooth structure. Failure to block-out may make tray removal difficult, or cause extraction of natural teeth or prosthesis.

**Impression Tray:**

All impression trays generally used for precision impressions.

- ▶ For sufficient adhesion, apply a thin layer of Polyether Adhesive to the tray and allow to dry completely according to the Polyether Adhesive Instructions for Use.

**Retraction**

Suitable retraction agents are e.g., aluminum hydrochloride- and aluminum sulfate-based solutions. Retraction with threads or rings containing epinephrine (adrenaline), 8-hydroxyquinoline sulphate or iron (III) sulphate may impair the setting of polyether impression materials.

- ▶ Keep the area dry from which an impression is to be taken.
- ▶ Hemostatic threads or rings may be used for subgingival preparations.
- ▶ Thoroughly remove the residue of the retraction agent by rinsing and drying.

**Dosing and Mixing**

- ▶ Dispense equal lengths of base paste and catalyst to a mixing pad. Over- and underdosing of the catalyst does not affect processing time, but impairs the quality of the impression.

- ▶ Mix pastes for 45 sec with a spatula. Touch off at the edge of the pad unmixed quantities of paste adhering to the spatula, and repeatedly take up and streak out until the color of the mixture is uniform. Do not mix by stirring in order to avoid the formation of bubbles!

**Times**

Mixing	Processing time from start of mixing	Setting time from start of mixing	Residence time in the mouth
min:sec	min:sec	min:sec	min:sec
00:45	02:45	06:00	03:15

These are the processing times at 74° F. For higher temperatures shorten, and lower temperatures prolong total processing times.

**Impression Taking**

**Mono-phase technique**

- ▶ By means of the filling device, place the mixture of pastes in the Elastomer syringe.
- ▶ Then load the tray previously prepared with adhesive.
- ▶ Apply the mixture into the sulcus or cavity from bottom up, always keeping the tip of the intra-oral tip immersed in the material and in contact with the surface of the tooth.
- ▶ Immediately after application, position the loaded tray in the mouth, and hold steady until setting is complete without applying pressure.
- ▶ To remedy initial adhesion (“setting the valve”), especially with upper jaw impressions, remove the tray from the gingiva on one side in a posterior position. If this proves difficult it may be necessary to carefully blow some air or water between the impression and the gingiva.

**Functional impression**

- ▶ After preparing the custom tray with adhesive, load the tray with paste mix, place it in the desired location in the mouth, and ask the patient to perform functional excursions.
- ▶ If required, use Permadyne Garant 2:1 or Permadyne low consistency for precise impression taking.

**After Impression Taking**

- ▶ Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

**Hygiene**

- ▶ Place the impression in a standard disinfectant solution commonly used for impressions for the period of time recommended by the manufacturer. Excessive disinfection may damage the impression.
- ▶ After disinfection, rinse the impression under running water for approx. 15 sec.

**Model Preparation**

- ▶ Prepare a cast from the impression with a commercial special stone plaster not earlier than 30 min and not later than 14 days after impression taking.
- ▶ To avoid introducing bubbles into the model, briefly rinse with water and dry with air prior to pouring. Do not use surfactants as these impair the quality of polyether impressions, and are not required!

- ▶ Polyether impressions can be silver-coated, whereas copper-coating is not feasible.

**Cleaning**

- ▶ Paste that has not set may be removed with ethanol or by rinsing with water and soap.
- ▶ The adhesive can be removed from metal trays with acetone.

**Notes**

- At temperatures below 18° C/64° F, the pastes cannot be extruded from the tube due to their increased viscosity. However, at room temperature processability is re-established without compromising quality.
- Direct exposure to sun and moist storage conditions damage the impression.
- Polyether impressions should not be exposed to solvent-containing liquids, as this may result in swelling and imprecise modelling.
- 3M ESPE Polyether materials may only be combined with other 3M ESPE polyether materials. Doing so the working and setting times of the individual products must necessarily be kept.

**Storage and Stability**

Store the product at 18–25° C/64–77° F. **Do not refrigerate!**

Do not use after the expiration date.

Store impressions dry and below 30° C/86° F in the dark.

**Customer Information**

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

**Warranty**

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH’s sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

**Limitation of Liability**

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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