GC Fuji I
Prior to use, carefully read the instructions for use.
EN

RADIOPAQUE GLASS IONOMER LUTING CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of metal-based inlays, onlays, crowns, bridges and posts.
- Cementation of high strength (zirconia based) all ceramic crowns and bridges.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder/liquid ratio (g/g)	1,8/1,0
Mixing time (min., sec.)	20"
Working time (min., sec.) (at 23°C) (from start of mixing)	2’00"
Final finishing commencing time (after seating the restoration)	4’30"

- TOOTH PREPARATION
 - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
 - Clean the prepared teeth with pumice and water.
 - Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).
- RESTORATION PREPARATION

Make sure that the restoration is pretreated and handled according to the manufacturer’s instructions.
- POWDER AND LIQUID DISPENSING
 - The standard powder to liquid ratio is 1.8g/1.0g (1 level scoop of powder to 2 drops of liquid).
 - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.
- MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, add all the powder to the liquid. Mix rapidly for 20 seconds. Note :

When mixing larger amounts, divide the powder into two equal parts. Mix the first portion with all of liquid for 5 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15 seconds (total 20 seconds).
- CEMENTATION
 - Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperatures shorten working time.
 - Maintain moderate pressure.
 - Start removing excess cement when the excess cement feels rubbery.
 - Finishing can be started 4 minutes 30 seconds after seating the restoration.

STORAGE

Store in a cool and dark place (4-25°C) (39.2-77.0°F). (Shelf life : 3 years from date of manufacture)

SHADE

Light yellow

PACKAGES

- 1-1 package : 35g powder, 25g (20mL) liquid, powder scoop, mixing pads (No. 20).
- Bottle of 35g powder with scoop.
- Bottle of 25g (20mL) liquid.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with components of other glass ionomer cements.
- This product is not indicated for filling or core build-up.

MANUFACTURED BY GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V. Research-Park, Haasroode-Leuven 1240, Interleuventaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC. 3737 West 127th Street, Alsip, IL 60803 U.S.A. TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD. 11 Tampines Concourse, #03-05, Singapore 528729 TEL: +65 6546 7588

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