



200/BOX



200/BOX

Perform Like A Superstar. Transform Your Game.

Treat your hands like superstars with Aurelia's newest comfort and performance series nitrile gloves. Introducing **Perform**™ and **Transform**™ — our newest nitrile gloves that combine all the best benefits of a nitrile glove in a proprietary formulation that take strength, comfort and performance to an all new level.

Perform™



- Our thinnest ever nitrile glove for ultimate comfort
- Greater tactile sensitivity
- Designed for optimal performance
- New teal color

Transform™




- Ultra thin for superior tactile sensitivity
- Combines long-lasting comfort, durability, and decreased hand fatigue
- Improved barrier protection
- New TransBlue™ color

For an experience beyond the ordinary,
visit aureliagloves.com and request a **FREE** sample kit.

1-877-AURELIA (1-877-287-3542) www.aureliagloves.com



Join us on Facebook at www.facebook.com/aureliagloves



AURELIA™
Sign of Comfort™

Treat your hands like superstars with Aurelia's newest comfort and performance series nitrile gloves. Introducing **Perform™** and **Transform™** — our newest nitrile gloves that combine all the best benefits of a nitrile glove in a proprietary formulation that take strength, comfort and performance to an all new level.



Aurelia Perform™

Length (mm/in): 240/9.5
 Cuff Thickness (mils): 2
 Palm Thickness (mils): 2.6
 Finger Thickness (mils): 3.6
 Tensile Strength:
 Before Aging (MPa): min 14
 After Aging (MPa): min 14
 Elongation:
 Before Aging (%): min 500
 After Aging (%): min 400
 Material: Nitrile (Latex Free)
 Glove Design : Ambidextrous
 Glove Feature : Powder Free
 and Finger Textured
 Color: Teal Green
 Packaging: 2000 gloves/case
 10 boxes 200 gloves/case

Item Numbers: XS-92395
 S-92396
 M-92397
 L-92398
 XL-92399

Aurelia Transform™

Length (mm/in): 240/9.5
 Cuff Thickness (mils): 2
 Palm Thickness (mils): 2.8
 Finger Thickness (mils): 4
 Tensile Strength:
 Before Aging (MPa): min 14
 After Aging (MPa): min 14
 Elongation:
 Before Aging (%): min 500
 After Aging (%): min 400
 Material: Nitrile (Latex Free)
 Glove Design : Ambidextrous
 Glove Feature : Powder Free
 and Finger Textured
 Color: Trans Blue
 Packaging: 2000 gloves/case
 10 boxes 200 gloves/case

Item Numbers: XS-98895
 S-98896
 M-98897
 L-98898
 XL-98899

1-877-AURELIA (1-877-287-3542) www.aureliagloves.com

 Join us on Facebook at www.facebook.com/aureliagloves



AURELIA
 Sign of Comfort™