

## SAFETY DATA SHEET

**1. Identification of the substance or mixture and of the supplier***1.1 Product Name:***OneGloss***1.2 Company /Undertaking identification*Manufacturer's Name: **SHOFU INC.**

Address: 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983 ,JAPAN

Phone: +81-75-561-1112

Fax: +81-75-561-2273

Section in Charge: Quality Assurance Section

Date Prepared: 2014.9.10

*1.3 Emergency Telephone Number*

+81-75-561-1112

**2. Hazardous identification**

## GHS CLASSIFICATION

It is not possible to classify it.

## LABEL ELEMENTS

SYMBOL	—
SIGNAL WORD	—
HAZARD STATEMENTS	—

## PRECAUTIONARY STATEMENTS

[Prevention]

Do not handle until all safety precautions have been read and understood.

[Response]

If you feel un-well, get medical advice/attention.

[Storage]

Store in a dry place. Take measures to avoid mechanical and chemical damage.

[Disposal]

Dispose of contents and container in accordance with regulation.

**3. Composition/information on ingredients**

Working part :	Synthetic rubber (Polyvinylsiloxane)
	Abrasive grain (Al <sub>2</sub> O <sub>3</sub> )
	Silicon oxide (SiO <sub>2</sub> )
	Others

No hazardous material is present in this grinding stone.

**4. First-aid measures**

- 4.1 Eye contact :* Wash immediately plenty of water and seek medical advice, since the grinding dust cause scratched harm to the cornea if contact.
- 4.2 Skin contact :* N/A
- 4.3 Ingestion :* N/A
- 4.4 Inhalation:* Remove from exposure .(if inhaled the dust from grinding)

## SAFETY DATA SHEET

**5. Fire-fighting measures**

- 5.1 Suitable Extinguishing Media:* CO<sub>2</sub>, Dry chemical, Foam
- 5.2 Extinguishing Media Should not be used :* N/A
- 5.4 Special Protective Equipment:* Wear fire protective cloth, safety goggles and dust mask.

**6. Accidental release measures**

- 6.1 Personal Precautions:* N/A
- 6.2 Methods for Cleaning Up:* Pick up and collect in a suitable container.

**7. Handling and storage**

- 7.1 Handling:* Use in accordance with instruction Regarding rotating speed.
- 7.2 Storage:* Store in a dry place. Take measures to avoid mechanical and chemical damage.

**8. Exposure controls/personal protection**

- | <i>8.1 Exposure limits:</i>        | Chemical name                         | CAS No. | ACGIH TLV |
|------------------------------------|---------------------------------------|---------|-----------|
|                                    | N/A                                   |         |           |
| <i>8.2 Respiratory Protection:</i> | Dust respirator and/or local exhaust. |         |           |
| <i>8.3 Hand Protection:</i>        | N/A                                   |         |           |
| <i>8.4 Eye Protection:</i>         | Safety goggles                        |         |           |
| <i>8.5 Skin Protection:</i>        | General protective clothing           |         |           |

**9. Physical and chemical properties**

- 9.1 Appearance/Odor :* Grayish-white colored silicone rubber bonded polishing stone (or points).
- 9.2 Boiling Point/Boiling Range:* N/D
- 9.3 Melting Point/Melting Range:* N/D
- 9.4 Flash point:* N/D
- 9.5 Auto flammability:* N/D
- 9.6 Explosive Properties:* N/D
- 9.7 Oxidizing Properties:* N/D
- 9.8 Vapour Pressure:* N/D
- 9.9 Relative Density:* N/D
- 9.10 Solubility:* *water solubility* Non-soluble

**10. Stability and reactivity**

- Stable in a normal temperature.
- 10.1 Condition to Avoid:* Avoid direct sunlight and high temperature.
- 10.2 Materials to Avoid:* N/A
- 10.3 Hazardous Decomposition Products:* Will not occur

**11. Toxicological information**

No information

**12. Ecological information**

No information

## SAFETY DATA SHEET

---

**13. Disposal considerations**

---

Dispose of contents/container to in accordance with local/regional/national/international regulations (to be specified).

---

**14. Transport information**

---

Airline restriction information ; It doesn't correspond.

UN No. ; It doesn't correspond.