

# Hygiene Heaven!

**Mobility. Choice. Control.**  
*YES, you can have it all!*

## new! AeroPro™

### The New Standard in Cordless Prophy Systems

Now there's finally a light-weight, ergonomic, cordless and pedal-less handpiece that gives you the mobility you want, the choice of prophy angles you like, and the control you need.

Based on extensive input from dental professionals, AeroPro's innovative design helps reduce hand fatigue by allowing you to rotate the angle 360° for improved intraoral access and comfort. Its quick-touch button enables fingertip control for easy switching between low and high speeds. Its quiet yet powerful motor maintains a consistent speed during polishing – no extra pressure required. ChargeSMART™ Battery Technology gives you long-lasting power to get through the workday. And yes, AeroPro™ is CDC-compliant with a convenient infection prevention protocol. Best of all, your patients get a great polish. Every time.



**premier** Inspired solutions  
for daily dentistry®

## new! AeroPro™ Cordless Prophy System

### Mobility

- ✓ Total mobility with no cord drag or rheostat
- ✓ ChargeSMART™ Battery Technology gets you through the day

### Choice

- ✓ Works with Premier 2pro® and most disposable prophy angles
- ✓ Quick-touch button to easily switch between speeds

### Control

- ✓ Ergonomic, lightweight, well-balanced design
- ✓ Rotatable prophy angle for 360° intraoral access and comfort
- ✓ Compliant to CDC infection prevention guidelines

**5500500 AeroPro™ Basic Kit**  
(1) Autoclavable Outer Sheath,  
(1) Bag of 48 2pro® Disposable  
Prophy Angles (Short/Soft),  
1 each: Motor Component, Charging Station,  
AC Adapter, Stand, Box of 500 Disposable Barriers,  
Accessories

**5500510 AeroPro™ Complete Package**  
(3) Autoclavable Outer Sheaths,  
(2) Bags of 144 2pro® Disposable  
Prophy Angles (Short/Soft),  
1 each: Motor Component, Charging Station,  
AC Adapter, Stand, Box of 500 Disposable Barriers,  
Accessories

*Best  
value!*

**5500520 AeroPro™ Autoclavable Outer Sheath**  
**5500530 AeroPro™ Disposable Barriers** - Box of 500



**Accepts most  
disposable angles**  
Can be rotated 360°

**Quick-touch button**  
Easily controls speeds

**Dimpled silicone grip**  
Improves tactile feel while  
preventing slippage

**Quiet and powerful  
German-made motor**  
Constant speed without  
extra pressure

**ChargeSMART™  
Battery Technology**  
Charges quickly and provides  
lasting power for all day use

**premier®** Inspired solutions  
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com