

JANUARY 2013

Dental

Your Source for Product Evaluations and Information

PRODUCT SHOPPER®

www.dentalproductshopper.com

OVER 200 PRODUCTS IN THIS ISSUE

TEMPORIZATION

DENMAT HOLDINGS, LLC **Snap-On Smile**



A practical and comfortable provisional solution that provides clinicians and patients with the esthetics they demand. **P.34**

We put the **SPOTLIGHT** on:

GLASS IONOMERS
SILMET
ProGlass Plus P.40



IMPRESSION TRAYS
CLINICAN'S CHOICE
Quad-Tray Xtreme P.52



MEET THE EVALUATORS

Dr. Denise Habjan, **P.64**



5 PEER-TO-PEER **EVALUATIONS:**

BENCO DENTAL
natural extensions
nitrile plus Fuchsia
P.68



Dental PRODUCT SHOPPER
BEST PRODUCT 2013



SILMET ProGlass Plus

Silmet's ProGlass Plus is a resin-modified glass ionomer luting cement. This glass ionomer features high bond strength and is designed to be used for the cementation of metal, ceramic inlays and crowns, inlays and onlays, porcelain-fused-to-metal and resin crowns, and bridges.

Designed with very low film thickness and excellent handling properties, ProGlass Plus can also change colors to assist in placement. It then shifts to a translucent A2 shade when set. Additionally, no etch is required before placement.

According to Silmet, ProGlass Plus does not cause any postoperative sensitivity and delivers strong adhesion to teeth and metal. These properties ensure both patient and clinician

satisfaction. The Bisphenol-free formulation of ProGlass Plus also delivers high radiopacity and fluoride release.

ProGlass Plus is offered in powder/liquid kits as well as in Silmet's patent-pending Luer Lock capsule with screwable tip. For ease of capsule insertion into the amalgamator, the application tip can be easily screwed onto the capsule after trituration, allowing easy and accurate application of the material to the designated surface.



SILMET LTD
FREE Product Information
CALL 800.697.7702 ext. 78200
CLICK www.dps.li/a/3AO-200
CIRCLE Reader Service 200