



Hemostatic Solution

SAFETY DATA SHEET

SECTION I - PRODUCT INFORMATION

PRODUCT NAME: Hemostat 30cc HEMOSTATIC SOLUTION

ITEM NO.: 25-03677

MANUFACTURER: Cargus International Inc.

135 North Route 9W

Congers, NY 10920

FOR INFORMATION CALL: (845) 267-2600

SECTION II - HAZARDOUS INGREDIENTS/IDENTITY INFORMATION

| <u>HAZARDOUS COMPONENTS</u> | <u>CAS NO.</u> | <u>OSHA PEL</u> | <u>ACGIH TLV</u> | <u>% (OPTIONAL)</u> |
|-----------------------------|----------------|-----------------|------------------|---------------------|
| ALUMINUM CHLORIDE | 7446-70-0 | 3.10g/Kg | 3.10g/Kg | 21% |
| PROPYLENE GLYCOL | 57-55-6 | 21g/Kg | 21g/Kg | - |

SECTION III - PHYSICAL DATA

BOILING POINT: 210 Degrees Farenheit

VAPOR PRESSURE (mm Hg): N/I

VAPOR DENSITY (AIR=1): N/I

SOLUBILITY IN WATER: Soluble

APPEARANCE AND ODOR: Clear, pale liquid, characteristic odor

SPECIFIC GRAVITY: 1.1

MELTING POINT: N/A

EVAPORATION RATE: N/I

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

FLASH POINT: N/A

FLAMMABLE LIMITS: Non-flammable **LEL:** N/A **UEL:** N/A

EXTINGUISHING MEDIA: Carbon Dioxide, water spray, foam or dry chemical

SPECIAL FIRE FIGHTING PROCEDURES: None

UNUSUAL FIRE AND EXPLOSION HAZARDS: None

SECTION V - REACTIVITY DATA

STABILITY: STABLE

CONDITIONS TO AVOID: None under anticipated normal use.

INCOMPATIBILITY (MATERIALS TO AVOID): None

HAZARDOUS DECOMPOSITION OR BYPRODUCTS: None

HAZARDOUS POLYMERIZATION: Will not occur.



SECTION VI – HEALTH HAZARD DATA

ROUTES OF ENTRY:

- **INHALATION:** YES - **SKIN:** YES - **INGESTION:** YES

HEALTH HAZARD (ACUTE AND CHRONIC): Mild skin irritation over prolonged exposure to skin, irritating to eyes.

CARCINOGENICITY: NTP: No IARC MONOGRAPHS: No OSHA REGULATED: No

SIGNS AND SYMPTOMS OF EXPOSURE: Mild Irritation.

MEDICAL CONDITIONS GENERALLY AGGRAVATED BY EXPOSURE: None Known.

EMERGENCY AND FIRST AID PROCEDURES:

- **EYES:** Flush with water for at least 15 minutes. Seek medical attention.
 - **SKIN:** Wash with soap and water. If irritation occurs, seek medical attention.
 - **INGESTION:** Consult a physician.
-

SECTION VII – PRECAUTIONS FOR SAFE HANDLING AND USE

STEPS TO BE TAKEN IN CASE MATERIAL IS RELEASED OR SPILLED: Wash with plenty of water, mop.

WASTE DISPOSAL METHOD: Follow applicable Federal, State and Local regulations.

PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE: Keep container tightly closed. Wash hands with soap and water after handling.

OTHER PRECAUTIONS: None under normal use.

SECTION VIII – CONTROL MEASURES

RESPIRATORY PROTECTION: None Required

VENTILATION - LOCAL EXHAUST: N/A

MECHANICAL: N/A

SPECIAL: N/A

OTHER: N/A

PROTECTIVE GLOVES: Recommended

EYE PROTECTION: Recommended

OTHER PROTECTIVE CLOTHING OR EQUIPMENT: None Required

WORK/HYGIENE PRACTICES: Wash hands with soap and water after each handling.

INFORMATION HEREIN IS GIVEN IN GOOD FAITH AS AUTHORITATIVE AND VALID. HOWEVER, NO WARRANTY, EXPRESSED OR IMPLIED CAN BE MADE, AND DEEPAK PRODUCTS INC. ASSUMES NO LIABILITY RESULTING FROM ITS USE.
