

# THE UNIVERSAL COMPOSITE THAT KEEPS POLISHING IN THE MOUTH



## G-ænial Sculpt®

Light-Cured Universal Nano-Hybrid Compactable Composite

- Excellent polish and high gloss retention
- Self-polishing effect enables additional polish in mouth
- Extraordinary wear resistance due to uniform nano-filler dispersion
- Easy to shape, yet stays in place
- High radiopacity (>300% Al)
- Pairs perfectly with G-Premio BOND™, our universal, 8th generation bonding agent






The unique surface treatment of nano-inorganic filler and high density and uniform dispersion technology. The technology achieves wear resistance and good physical properties (shining surface texture like natural enamel with lasting esthetics). Even in areas that are difficult to polish, filled G-ænial Sculpt® obtains gloss with daily brushing.



Courtesy of Dr. Robert Marus.



G-ænial Sculpt®'s fine, 300-nanometer barium glass particle allows for high strength esthetics, excellent polish and high gloss retention. Its low shrinkage stress and ideal handling properties allow for placement in all classes of restorations. The high filler loading (78% wt) and homogenous dispersion of nano-filler maximizes polish retention and wears better than other universal composite restoratives.

### G-Premio BOND™ is a Great Companion to G-ænial Sculpt®

G-Premio BOND™ is a one-component, truly universal adhesive bonding agent, which is flexible for use in total, selective and self-etch techniques. Offering excellent stability, G-Premio BOND™ has an extended working time of up to 5 minutes. The product bottle is designed for more precise dispensing control and features a silicone tip to help eliminate waste. G-Premio BOND™ Dual Cure Activator (DCA) provides additional flexibility by accommodating a self-cure option. G-Premio BOND™ may also be used as a desensitizer.



- Flexible for use in total, selective and self-etch bonding techniques
- Fast application (10 seconds)
- Good stability with extended working time (up to 5 minutes)
- New bottle design for better dispensing control

#### G-ænial Sculpt® Syringe Refills

Contains: One Syringe (2.0mL).

009152 A1	009153 A2	009154 A3	009155 A3.5	009156 A4
009168 AE	009161 AO1	009162 AO2	009163 AO3	009157 B1
009158 B2	009159 B3	009164 BW	009160 C3	009166 CV
009167 JE	009165 XBW			

#### G-ænial Sculpt® 20-Unit Refills

Contains: 20 unitips (0.16 mL each).

009169 A1	009170 A2	009171 A3	009172 A3.5	009173 A4
009174 B1	009175 B2	009176 B3	009177 C3	

#### G-ænial Sculpt® 10-Unit Refills

Contains: 10 unitips (0.16 mL each).

009185 AE	009178 AO1	009179 AO2	009180 AO3	009180 BW
009183 CV	009184 JE	009182 XBW		

#### 012691 G-Premio BOND™ Kit

Contains: One 5mL bottle, 50 disposable applicators (fine), 20 disposable dispensing dishes.

#### 012696 G-Premio BOND™ 5mL Bottle Refill

#### 009275 G-Premio BOND™ 50 Unit Dose (0.1mL each) Refill

#### 009552 G-Premio BOND™ Dual Cure Activator (DCA) 3mL Bottle Refill