

3M Science.
Applied to Life.™

**Efficient Restorative
Procedure Solution—
Repair & Protect**

**Because they like
staying on schedule
as much as you do.**



Your time matters as much as theirs ...

When your schedule gets thrown off—as we all know can happen—patients can't hide their disappointment, and your day runs long. So, to help you stay on schedule, we've simplified the procedure you perform most frequently—the posterior restoration.

At the core of the procedure is 3M™ Filtek™ One Bulk Fill Restorative ... with increased opacity for greater esthetics. It can be placed in one fast, easy increment up to 5 mm—without the need for layering or additional expensive dispensing devices that can slow you down.



By using it with our other 3M products that are easy to use, you'll move through posterior restorations with speed and simplicity—saving time without compromising esthetic results.

Your patients will be happy—and so will you when you can end your day on time.

Efficient and consistent procedure for Class II restorations

- ▶ Simple, one-step bulk placement up to 5 mm with Filtek One Bulk Fill Restorative
- ▶ High bond strength and virtually no post-op sensitivity with 3M™ Scotchbond™ Universal Adhesive
- ▶ Deep, uniform cure with the 3M™ Elipar™ DeepCure-S LED Curing Light
- ▶ Polish restorations to a high-gloss finish with the 3M™ Sof-Lex™ Diamond Polishing System
- ▶ Protect with 3M™ Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste



Repair and protect

Fast and easy from start to finish.




3M™ Scotchbond™ Universal Adhesive

Versatile. High bond strength. One simple technique and quick 35-second application. Virtually no post-operative sensitivity.



3M™ Filtek™ One Bulk Fill Restorative

Simplicity of bulk placement up to 5 mm deep. Low wear for durability and stress relief. Optimized opacity for greater esthetics.



3M™ Elipar™ DeepCure-S LED Curing Light

Predictable, reliable cure. Optimized light guide is more comfortable for patients and provides easier access to posterior teeth.



3M™ Sof-Lex™ Diamond Polishing System

Easily adapts to irregular, convex and concave tooth surfaces and imparts paste-like gloss with the convenience of a rubberized system. Fast and easy to use.



3M™ Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste

3M's Tri-Calcium Phosphate (TCP) helps repair white spots and protect against cavities. It replaces the toothpaste in the normal brushing routine, making compliance easy.



Prescription only. Use in adults and children 6 years of age and older.

Use once daily in place of conventional toothpaste. Go to www.3M.com/dental to obtain the full prescribing information or call 1-800-634-2249. You are encouraged to report negative side effects of prescription drugs to the FDA. Visit www.fda.gov/medwatch, or call 1-800-FDA-1088.

Ordering Information



3M™ Filtek™ One Bulk Fill Restorative

Refill SHADE	Item # SYRINGE	Item # CAPSULE
A1	4866A1	4867A1
A2	4866A2	4867A2
A3	4866A3	4867A3
B1	4866B1	4867B1
C2	4866C2	4867C2

SYRINGE SHADE REFILLS include: 1 – 4g syringe; Technique Guide; Instructions for Use.

CAPSULE SHADE REFILLS include: 20 – 0.2g capsules; Technique Guide; Instructions for Use.



3M™ Scotchbond™ Universal Adhesive

Item # Product Information

41254	3M™ Scotchbond™ Universal Adhesive Intro Kit—Vial Kit includes: 1 – 5ml adhesive vial; 50 disposable mixing wells; 50 disposable microbrushes; 1 – 3ml etchant syringe
41258	3M™ Scotchbond™ Universal Adhesive Refill Includes: 1 – 5ml adhesive vial
41255	3M™ Scotchbond™ Universal Adhesive Intro Kit—Unit Dose 40 Includes: 40 unit dose; 1 – 3ml etchant syringe
41256	3M™ Scotchbond™ Universal Adhesive Intro Kit—Unit Dose 100 Includes: 100 unit dose; 1 – 3ml etchant syringe
41257	3M™ Scotchbond™ Universal Adhesive Intro Kit—Unit Dose 200 Includes: 200 unit dose; 1 – 3ml etchant syringe
41263	3M™ Scotchbond™ Universal Etchant Etching Gel Includes: 2 – 3ml etchant syringes; 50 dispensing tips
41264	3M™ Scotchbond™ Universal DCA Dual Cure Activator Includes: 1 – 5ml vial



3M™ Elipar™ DeepCure-S LED Curing Light

Item # Product Information

76976	3M™ Elipar™ DeepCure-S LED Curing Light Includes: Handpiece (Cordless), Charging Base (230V); Li-ion Battery; 10mm Light Guide; Eye Shield
76981	3M™ Elipar™ DeepCure-S Light Guide, 10mm
76984	3M™ Elipar™ DeepCure Eye Shield
76985	3M™ Elipar™ DeepCure-S Rechargeable Li-ion Battery

3M™ Sof-Lex™ Diamond Polishing System

Item # Product Information



5092-IM	3M™ Sof-Lex™ Diamond Polishing System—Introductory Kit Kit includes: 5 – 3M™ Sof-Lex™ Pre-Polishing Spirals (Beige); 5 – 3M™ Sof-Lex™ Diamond Polishing Spirals (Pink); 1 – 3M™ Sof-Lex™ RA Mandrel; Instructions for Use; Technique Card
5090	3M™ Sof-Lex™ Pre-Polishing Spirals—Refill Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige)
5091	3M™ Sof-Lex™ Diamond Polishing Spirals—Refill Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink)



3M™ Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste

Item # Product Information

12115	Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste—Vanilla Mint Flavor – 4oz tube
12115BG	Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste—Bubble Gum Flavor – 4oz tube
12115SM	Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste—Spearmint Flavor – 4oz tube



www.3M.com/EfficientRestorations

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M, ESPE, Clinpro, Elipar, Filtek, Scotchbond and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA. © 3M 2016. All rights reserved.
70-2013-0805-6