



Polycarbonate Crowns

FEATURES & BENEFITS

- Contour and crimp similar to metal crowns for easy adaptation to margins
- Universal shade can be adjusted using cements and liners
- Can be used for short or longer-term use providing the option of offering an economical aesthetic long-term temporary
- Universal shade and numbering
- Will not crack, split, shatter or chip
- Easily scissor trimmed, contoured, crimped and polished
- Thin design prevents gingival irritation
- May be cemented with any plastic, zinc oxide or orthophosphate cement

UPPER CENTRAL Right

Size	Item#
100	P100
101	P101
10	P10
11	P11
12	P12
13	P13
14	P14

UPPER CENTRAL Left

Size	Item#
15	P15
16	P16
17	P17
18	P18
19	P19
102	P102
103	P103

RIGHT CUSPID Upper & Lower

Size	Item#
300	P300
301	P301
30	P30
31	P31
32	P32
33	P33
34	P34

LEFT CUSPID Upper & Lower

Size	Item#
35	P35
36	P36
37	P37
38	P38
39	P38
302	P302
303	P303

UPPER LATERAL Right

Size	Item#
2	P2
20	P20
21	P21
22	P22
23	P23
24	P24

UPPER LATERAL Left

Size	Item#
25	P25
26	P26
27	P27
28	P28
29	P29
200	P200

LOWER ANTERIOR Long

Size	Item#
60	P60
61	P61
62	P62
63	P63
64	P64

LOWER ANTERIOR Short

Size	Item#
69	P69
68	P68
67	P67
66	P66
65	P65

1st BICUSPID

Size	Item#
40	P40
41	P41
42	P42
43	P43
44	P44

2nd BICUSPID

Size	Item#
54	P54
53	P53
52	P52
51	P51
50	P50

C180.....Crown Kit 180 cs - 3 of each size
Replacement Crowns 5 pcs per box