

Silmet's Selection of Glass Ionomer Cements

Features:

- Bisphenol FREE formulation
- High Radiopacity
- Excellent handling properties
- No etch required prior to placement
- Good biocompatibility
- Fluoride release

Glass Ionomer Capsule

Tailored to offer a wide range of restoration and cementation solutions, Silmet's new Innovative Luer Lock Cap GI Capsule, features a screwable tip, to increase the ease of use and achieve best results during the use of ProGlass GIC line of products. The Patent pending Luer Lock cap enables to easily remove and attach the capsule tip for easy insertion into amalgamator.



SILMET
www.silmetdental.com

For more information,
call our Support Center:
1-800-972-7917/866-300-3657

	ProGlass One	ProGlass Plus	ProGlass Nine	ProGlass Two	ProGlass Two LC	ProGlass Silver	ProBase
Powder/Liquid Ratio (g/g)	2.4/ 1.0 1 Red spoon 3drops	1.8 /1.0 1 Orange spoon 4 drops	4.1 / 1.0 1 Blue spoon 2 drops	3 / 1.0 1 Red spoon 2 drops	3 / 1.0 1 Blue spoon 2 drops	4.1 / 1.0 1Red spoon 1 drops	1.4 / 1.0 1 Orange spoon 4 drops
Working Time (min.sec)	2'30" – 3'	3'	2'30"	3'	3' Light Cure for 20" (1.8mm)	1'50"	3' Light Cure for 30" (1.8mm)
Setting Time (min.sec)	3'10"	3'	3'30"	3'50"	3'	4'	3'
Compressive Strength 24 hr. (MPa)	200	220	260	165	260	135	220
Diametrical Tensile strength 24hr. (MPa)	30.5	35	30	33	40.5	25.5	35

Restorative Type II / R

ProGlass Nine

Radiopaque Packable Glass Ionomer Cement

Similar to:
GC Fuji IX

ProGlass Two / Two LC

Radiopaque Glass-ionomer Restorative cement in two variations: Self Cure or Light Cure

Similar to:
Fuji II (GC)
Vidrion R (SS White)

Fuji IILC (GC)
Ketac Nano (3M ESPE)

Luting / Cementation Type I / C

ProGlass One

Glass Ionomer Luting Cement

Similar to:
Ketac Cem (3M ESPE)
Fuji I (GC)
Meron AC (VOCO)

ProGlass Plus

Resin Reinforced Glass Ionomer cement

Similar to:
Fuji PLUS (GC)
Fuji CEM (GC)
Rely-X / Vitremer (3M ESPE)

Core Build Up Type VII / N

ProGlass Silver

Silver Reinforced Glass Ionomer Cement for Core Build Up

Similar to:
Miracle Mix (GC)
Ketac Silver (3M ESPE)

Liner / Base Type III / F

ProBase

Light Cure, Glass Ionomer Liner/Base

Similar to:
Fuji LINING™ LC (GC)
Vitrebond (3M ESPE)

Presentation:

Powder/Liquid Working Kit
Trial Size Powder/Liquid Kit (approx. 1/2 size)
GI Capsules (Box of 50)

