



Lock-Tight Plastic Air/Water Syringe Tips with Rubber Seal



Water Channel



Rubber O-Ring Seal



MARK3 Seal Tight Replacement Air Water Syringe Tips are disposable tips with an inner yellow O-ring seal which provides dry air with every use.

FEATURES & BENEFITS

- Easy to remove and replace
- Eliminates risk of cross contamination
- Eliminates sterilization time
- Reduces risk of moisture damage during procedures
- Flexible, can bend to a 90° angle without affecting water flow
- Rigid, retracts lips and cheek without bending

The MARK3 Line
was created to provide
Premium Products
at an Affordable Price

White Tips

#1205 200/bag

#1206 1500/bag

Multicolored Tips

#1207 200/bag

#1208 1500/bag

Compare to Seal Tight by Kerr



Authorized for sale in Canada