

MARK3TM
THE SMART ALTERNATIVE



Disposable

Dental Needles

INDICATIONS

Routine administration of local dental anesthetics.

FEATURES & BENEFITS

- Less trauma to the tissue achieved through a triple-bevel needle, a polished & siliconized cannula and excellent flexibility of the needle.
- High grade surgical stainless steel tubing to reduce risk of breakage.
- Bevel mark for optimal injection position.
- Compares to Septodont, Sultan and other private label brands.

The MARK3 Line
was created to provide
High Quality Products
at an Affordable Price

order#	sizes
1630327 Gauge - 1" Short
1630427 Gauge - 1 ³ / ₈ " Long
1630530 Gauge - 3/ ₈ " Extra Short
1630630 Gauge - 1" Short

100 per box