



Make it easy on yourself!

Monobond Plus

**One Primer for
all restorative materials**

- Glass-ceramics
- Noble-metal alloys
- Zirconium-oxide ceramics
- Aluminium-oxide ceramics
- Base-metal alloys
- Composites (indirect restorations)



ivoclar
vivadent
passion vision innovation

Monobond Plus

One Primer

For the first time ever, a primer is available which you can use to condition **all types of restoration surfaces**: Monobond Plus.

The innovative combination of three different functional groups – silane methacrylate, phosphoric acid methacrylate and sulfide methacrylate – enables you to establish a **strong and durable bond** to any dental restorative material.

One reaction time

Monobond Plus requires the same reaction time on all the various material surfaces. This makes it even easier for you to cement restorations:

- Prepare the restoration surface
- Apply Monobond Plus and **allow it to react for 60 seconds**
- Disperse any excess with air

Done!

Easy storage

No special storage conditions are required, as Monobond Plus can be stored **at room temperature**.

Moreover, Monobond Plus helps you streamline your inventory: You only require **one primer** for all the different restorative materials.



Delivery form

Refill

1 bottle Monobond Plus

5 g



**Monobond Plus makes
your life a little easier –**

**AS FAR AS THE CEMENTATION OF
RESTORATIONS IS CONCERNED!**

Descriptions and data constitute no warranty of attributes.

Printed in Liechtenstein
© Ivoclar Vivadent AG, Schaan / Liechtenstein
627562/0309/e/BVD

Ivoclar Vivadent AG
Bendererstr. 2
FL-9494 Schaan
Fürstentum Liechtenstein
Tel. +423 / 235 35 35
Fax +423 / 235 33 60
www.ivoclarvivadent.com

ivoclar
vivadent[®]
passion vision innovation