

**MARK3**<sup>®</sup>  
THE SMART ALTERNATIVE



4.3

Dental **PRODUCT SHOPPER**

**BEST PRODUCT**

**MADE IN USA**



"I always thought the gold standard was Luxatemp, but this might unseat it as king."

*Dr. Tiffany Lee, DDS  
Torrance, CA*

Evaluator quote from  
Dental Product Shopper's  
Product Review, June 2021.

## 10:1 Automix

# Temporary Crown & Bridge

Self-Curing Temporary Crown & Bridge Material. Unique formulation features excellent mechanical and esthetic properties for short & long-term provisionals.

## FEATURES & BENEFITS

- Exceptional strength and abrasion resistance
- Defined elastic phase for easy removal
- True color stable Vita shades with natural fluorescence
- Minimal heat (= body temp.) during polymerization
- Easy to polish to a long lasting luster

## Technical Data

- Dosage and application time: Approx. 45 seconds
- Setting time in mouth: Approx. 2-3 minutes (hard elastic phase, from beginning of mixing)
- Setting time: Approx. 6-7 minutes (until setting time is completed)
- Mixing ratio: 10:1

1 Cartridge of 50m plus 10 Mixing Tips

**4101**.....Shade A1

**4102**.....Shade A2

**4103**.....Shade A3

**4105**.....Shade B1

Used with...

Dispensing  
Gun - 10:1

order# 110411  
Type 50



*Compare to Luxatemp*