



## Metal Core

# Quick Tips AIR/WATER SYRINGE TIPS

Quick Tips Disposable Metal Air Water Syringe Tips are Adec Style tips that eliminate the risk of cross-contamination, reduces sterilization time and risk of moisture damage during procedures.

### FEATURES & BENEFITS

- Interchangeable with most metal tips
- No conversion kits or adapters necessary
- Safe and easy to use
- Separate stainless steel water tube prevents air-water crossover
- Interchangeable with A-dec and DCI style syringes

**#1409** 100 tips per bag. 100 bags per case.

**#1411** 1600 tips per bag. 6 bags per case.



Authorized for sale in Canada

The MARK3 Line  
was created to provide  
Premium Products  
at an Affordable Price

*Compares to: Safe Tips EZ and Fast Tips*