

Surgical Gown Pack BCE-100



COMPOSITION	<ul style="list-style-type: none"> ✓ SMS light blue fabric 50 grams. ✓ Non-woven 100%. ✓ Disposable.
STERILIZATION	<ul style="list-style-type: none"> ✓ Gamma Radiation.
CHARACTERISTICS	<ul style="list-style-type: none"> ✓ Impervious to penetration of liquids and fluids, anti-glare, not transparent, anti-static and resistant to stress in normal use. ✓ Sterile and disposable. ✓ Adjustable cuffs.
CAUTIONS	<ul style="list-style-type: none"> ✓ For single use only. ✓ Sterilization not guaranteed if opened.
STORAGE	<ul style="list-style-type: none"> ✓ Keep in a cool dry place. ✓ Expiration date: 5 years since sterilization date.
PACKAGING	<ul style="list-style-type: none"> ✓ Pack containing: 1 Gown, 1 Hand towel, 1 Cover ✓ Carton with 50 gowns ✓ 16 boxes per pallet

Looking for **safety**, think of