

MARK3[®]
THE SMART ALTERNATIVE



Outstanding
for Successful
Tissue
Management



PRE-BENT Infusor Tips

FEATURES & BENEFITS

- Blunt, bent 19 gauge cannula with padded brush enables gentle pressure on the sulcus
- Universal luer lock style hub
- Single use disposable tip
- Curved cannula for optimal entry
- Soft bristle end allows for non-invasive brushing of solutions
- Places hemostatic agents precisely
- Effectively removes superficial coagulum

The MARK3 Line
was created to provide
High Quality Products
at an Affordable Price

#8025 100 per pack
100 packs per case

Compares to:
Stat-flo Tips by Vista Dental Products



AUTHORIZED FOR SALE IN CANADA