

# Supergel® Impression Material

Supergel® is a premium blended alginate material that mixes readily to a smooth, creamy consistency. Its slow dehydration and low surface tension provide consistent, accurate, bubble-free detailed reproductions every time. Pleasing sassafras flavor, controlled flow, heavy body, and gag reflex reduction provide greater patient acceptance.

Fluff powder by rotating or shaking bag before using. Use 72°F (23°C) water. If necessary to delay the pouring of the model, wrap a wet towel around the impression and store in a plastic bag or sealed container.



## ORDERING INFORMATION

### SUPERGEL®

- 0921825** Regular Set, Can - 1lb (454g)
- 0921826** Fast Set, Can - 1lb (454g)
- 0921818** Regular Set, Pouch - 1lb (454g)
- 0921819** Fast Set, Pouch - 1lb (454g)

## FEATURES AND BENEFITS

- 1 Smooth, creamy consistency
- 2 Bubble-free reproductions
- 3 Sassafras flavor
- 4 Gag reflex reduction