

SURFACE
DISINFECTANTS
BROCHURE



The most commonly used disinfectants are EPA registered wipes, ready-to-use spray and concentrated liquids.

There are four key features healthcare providers should consider when choosing a proper surface disinfectant:



CONTACT TIME



BROAD-SPECTRUM PATHOGEN CLAIMS



CLEANING PROPERTIES



MATERIAL COMPATIBILITY

The CaviCide™ and CaviWipes™ brands are designed for reducing cross-contamination in the healthcare community. Trust the surface disinfectant brand used by thousands of healthcare facilities everyday.



CaviWipes™ Bleach

CaviWipes Bleach is made for hard non-porous surfaces. These new surface disinfectant cleaning wipes are gentle on surfaces, yet tough on pathogens.

- Meets the CDC recommended 1:10 bleach dilution
- No Personal Protective Equipment (PPE) required¹
- Kills all product labeled organisms in 3 minutes including *C. diff*¹, TB, drug resistant bacteria, viruses and fungi.
- Compatible with most medical device equipment²
- One-step cleaner and disinfectant
- Low odor

Part #	Description	Packaging
13-9100	CaviWipes Bleach (6" x 10.5"), 90 count	12 canisters/case
13-1175	CaviWipes Bracket	12 brackets/case

Metrex offers inservicing and wall charts for your education needs.

References
 1. Special instructions for cleaning prior to disinfection against *Clostridium difficile* spores, HIV-1, HBV, HCV.
 2. See medical device equipment manual for compatibility information.
 3. When no gross debris is present.



CaviWipes1™ & CaviCide1™

CaviWipes1 and CaviCide1 are the only 1-minute low-alcohol disinfectants in the market. They are multi-purpose surface cleaners and disinfectants designed for every day use.

- 1-minute contact time for virucidal, bactericidal (including TB), and pathogenic fungi activity
- 1-step cleaner and disinfectant³
- Fragrance free & Bleach free
- Multi-purpose cleaner for every day use

Part #	Description	Packaging
13-5100	CaviWipes1 (6" x 6.75"), 160 count	12 canisters/case
13-5150	CaviWipes1 EXTRA LARGE (9" x 12"), 65 count	12 canisters/case
13-5155	CaviWipes1 EXTRA LARGE Singles (9" x 12"), 50 count	6 boxes/case
13-5224	CaviWipes1 Flat Pack (7" x 9"), 45 count	20 packs/case
13-1175	CaviWipes Bracket	12 brackets/case
13-5024	CaviCide1 - 24 ounce spray	12 bottles/case
13-5000	CaviCide1 - 1 gallon	4 bottles/case
13-5025	CaviCide1 - 2.5 gallon	2 bottles/case



CaviWipes™ & CaviCide™

CaviWipes and CaviCide are intermediate level disinfectants recommended for use on hard, non-porous surfaces in healthcare settings.

- Effective against TB and MRSA in 3 minutes
- Effective against HIV-1, HBV, and HCV in 2 minutes
- All in 1 cleaner and disinfectant
- Available in multiple configurations to meet your needs

Part #	Description	Packaging
13-1100	CaviWipes (6" x 6.75") - 160 wipes per canister	12 canisters/case
10-1090	CaviWipes (6" x 6.75") - 220 wipes per canister	12 canisters/case
13-1150	CaviWipesXL (9" x 12") - 65 wipes per canister	12 canisters/case
13-1155	CaviWipesXL (9" x 12") - 50 wipes per box	6 boxes/case
13-1224	CaviWipes in a Flat Pack (7" x 9") - 45 wipes per pack	20 packs/case
13-1175	CaviWipes Bracket	12 brackets/case
13-1024	CaviCide - 24 ounce spray	12 bottles/case
13-1000	CaviCide - 1 gallon	4 bottles/case
13-1025	CaviCide - 2.5 gallon	2 bottles/case

