

MARK3[®]
THE SMART ALTERNATIVE



Disposable - One Step

Finishing and Polishing System

The MARK3 Finishing & Polishing System is a single-use, 1-step polishing system designed to provide an outstanding polished surface on all direct-resin restoratives. It includes advanced nano composites.

FEATURES & BENEFITS

- Unique silicone rubber impregnated with aluminum oxide particles is 100% Latex Free
- Delivers perfect control to reduce surface roughness & bacterial retention
- Achieves superbly smooth finish and high gloss
- More durable & flexible than the leading brands
- Reinforced mandrel allows for optimal concentric torque
- Available in 3 shapes to enable access to interproximal and occlusal surfaces
- A very cost-efficient single-step polishing system



The MARK3 Line
was created to provide
High Quality Products
at an Affordable Price

#6340 Intro Kit - 10 x Point, 10x Disc, 10x Cup

#6341 - 30x Cup

#6342 - 30x Disc

#6343 - 30x Point



Compare to Pogo by Dentsply Caulk



AUTHORIZED FOR SALE IN CANADA