



480 South Democrat Road, Gibbstown, New Jersey, 08027 **1-800-333-3131**
Werner-von-Siemens Str. 14a D-78224 Singen, Germany **+49 7731 912101**



keystoneind.com 1.800.333.3131



TRAY SOLUTIONS



TRAY AWAYS®

Disposable Impression Trays

DESCRIPTION

Tray Aways® – The original plastic tray! Updated throughout the years, the trays can be used not only in dental, but also with implants.

Tray Aways® are the disposable impression trays of choice suited for use with alginate or other types of impression materials where detailed reproduction and dimensionally-stable results are vital. Made from rigid, high-impact FDA approved styrene, Tray Aways® assist in generating a consistent, unsurpassed impression free of distortion or bite through. Smooth, rounded edges, interior retention slots, and vertical perforation enhance material retention. The best value on the market.

SIZES

S M L


12 per bag

MATERIAL

Styrene

TYPE

Solid
Perforated



CONSISTENCY BITE REGISTRATION

Impression Material

DESCRIPTION

- Gives exceptional accuracy
- Plaster-like hardness
- Dependable, no-slump consistency

ITEM

- #6081000 Fast Set – 50ml Cartridge
- #6081100 Regular Set – 50ml Cartridge

*each includes (2) 50ml cartridges & mixing tips



CONSISTENCY VPS

Impression Material

DESCRIPTION

- Hydrophilic & strong
- Final impression material or alginate substitute
- All available viscosities are accurate

ITEM

- #6081200 Heavy Body Fast Set
- #6081300 Heavy Body Regular Set
- #6081400 Medium Body Fast Set
- #6081500 Medium Body Regular Set
- #6081600 Light Body Fast Set
- #6081700 Light Body Regular Set
- #6081900 HP Dispensing Gun, Split Cartridge
- #6081800 HP Mixing Tips, 48/Pkg

*each includes (4) 50ml cartridges



TACTM

Tray Adhesive Compound

DESCRIPTION

- Available in spray can or brush-on liquid
- Excellent for orthodontic use
- Securely holds alginate in place when impression is removed from the mouth

ITEM

- #0921880 Spray – 3.5 oz. (99 g)
- #0921881 Standard Kit (Brush-on Liquid)
- #0921882 Bulk Bottle (Brush-on Liquid) – 4 oz. (118 ml)





SUPERGEL® FRESH
Dustless Alginate Impression Material

DESCRIPTION

- Dust-controlled formula for less airborne particles when opening & mixing
- Blue in color; berry-cherry flavored
- Choice of setting speeds

ITEM

- #0921822** Regular Set, Can – 1lb (454g)
- #0921823** Fast Set, Can – 1lb (454g)
- #0921827** Regular Set, Pouch – 1lb (454g)
- #0921830** Fast Set, Pouch – 1lb (454g)

SUPERGEL®
Alginate Impression Material

DESCRIPTION

- Accurate, bubble-free reproductions
- Sassafras flavor
- Mixes to a smooth, creamy consistency
- Gag reflex reduction

ITEM

- #0921825** Regular Set, Can – 1lb (454g)
- #0921826** Fast Set, Can – 1lb (454g)
- #0921818** Regular Set, Pouch – 1lb (454g)
- #0921819** Fast Set, Pouch – 1lb (454g)



GREEN FLEXIBOLES
Mixing Bowls

DESCRIPTION

- Plaster, stone, investments & alginates will not stick
- Resistant to most chemicals, solvents, & oils
- **NOT** suitable for mixing acrylics
- Clean by flexing & wiping with a damp cloth

ITEM

- #1150060** Small
- #1150070** Medium
- #1150080** Large
- #1150090** X-Large

Perforated



#1 UPPER

SIZE
Large

ITEM
#0921884

(W) 3.06", 77.6 mm
(H) 2.41", 61.2 mm
(D) 0.88", 22.4 mm

#2 LOWER

SIZE
Large

ITEM
#0921887

(W) 3.02", 76.7 mm
(H) 2.30", 58.3 mm
(D) 0.83", 21.1 mm

#3 UPPER

SIZE
Medium

ITEM
#0921885

(W) 2.91", 73.9 mm
(H) 2.21", 56.2 mm
(D) 0.83", 21.3 mm

#4 LOWER

SIZE
Medium

ITEM
#0921888

(W) 2.93", 74.3 mm
(H) 2.18", 55.2 mm
(D) 0.79", 20.0 mm

#5 UPPER

SIZE
Small

ITEM
#0921886

(W) 2.66", 67.9 mm
(H) 2.00", 50.9 mm
(D) 0.85", 21.9 mm

#6 LOWER

SIZE
Small

ITEM
#0921889

(W) 2.78", 70.6 mm
(H) 2.06", 52.3 mm
(D) 0.77", 19.6 mm

Perforated Partial



7



8



9



10

#7 UR, LL

SIZE
Partial

ITEM
#0921894

(W) 2.24", 56.9 mm
(H) 2.29", 58.3 mm
(D) 0.79", 20.1 mm

#8 UL, LR

SIZE
Partial

ITEM
#0921895

(W) 2.15", 54.5 mm
(H) 2.33", 59.3 mm
(D) 0.86", 21.8 mm

#9 LOWER

SIZE
Anterior

ITEM
#0921896

(W) 2.56", 65.0 mm
(H) 1.50", 38.1 mm
(D) 0.75", 19.0 mm

#10 UPPER

SIZE
Anterior

ITEM
#0921897

(W) 2.52", 64.1 mm
(H) 1.47", 37.2 mm
(D) 0.76", 19.2 mm

Partial & Twist



ITEM	PARTIAL SOLID
#0921447	#30 Upper L & Lower R
#0921448	#31 Upper R & Lower L
#0921449	#32 Anterior Lower
#0921450	#33 Anterior Upper
#0921451	Set 110 contains 4 trays: #s 30, 31, 32, 33

ITEM	PARTIAL PERFORATED
#0921452	#30 Upper L & Lower R
#0921453	#30 Upper R & Lower L
#0921454	#32 Anterior Lower
#0921455	#33 Anterior Upper
#0921456	Set 109 contains 4 trays: #s 30, 31, 32, 33

ITEM	TWIST TRAY
#0921866	Twist Tray Perforated
#0921868	Twist Tray Solid

Perforated

ITEM	EDENTULOUS
#0921337	#1 Upper
#0921338	#2 Upper
#0921339	#3 Upper
#0921340	#4 Upper
#0921341	#5 Lower
#0921342	#6 Lower
#0921343	#7 Lower
#0921344	#8 Lower
#0921307	Set 103 contains 8 trays: 4U, 4L

ITEM	REGULAR
#0921654	#1 Upper
#0921668	#2 Upper
#0921475	#3 Upper
#0921476	#4 Upper
#0921663	#5 Upper
#0921664	#1 Lower
#0921665	#2 Lower
#0921477	#3 Lower
#0921478	#4 Lower
#0921666	#5 Lower
#0921479	#0 Upper (XL)
#0921180	#0 Lower (XL)
#0921667	Set 106 contains 10 trays: 5U, 5L

ITEM	PEDO/ORTHO
#0921872	#15 Upper
#0921873	#16 Upper
#0921874	#17 Upper
#0921875	#18 Upper
#0921876	#15 Lower
#0921877	#16 Lower
#0921878	#17 Lower
#0921879	#18 Lower
#0921718	Set 108 contains 8 trays: 4U, 4L



7s



8s



9s



10s

#7s UR, LL

SIZE
Partial

ITEM
#0921890

(W) 2.24", 56.9 mm
(H) 2.29", 58.3 mm
(D) 0.79", 20.1 mm

#8s UL, LR

SIZE
Partial

ITEM
#0921891

(W) 2.15", 54.5 mm
(H) 2.33", 59.3 mm
(D) 0.86", 21.8 mm

#9s LOWER

SIZE
Anterior

ITEM
#0921892

(W) 2.56", 65.0 mm
(H) 1.50", 38.1 mm
(D) 0.75", 19.0 mm

#10s UPPER

SIZE
Anterior

ITEM
#0921893

(W) 2.52", 64.1 mm
(H) 1.47", 37.2 mm
(D) 0.76", 19.2 mm

Solid Partial

EDENTULOUS TRAY AWAYS®

Disposable Impression Trays

DESCRIPTION

Edentulous Tray Aways® are single-use, perforated edentulous impression trays that can be easily drilled for implant impressions. These unique and economical trays allow for maximum material retention and accuracy without the use of an adhesive.

SIZES

S M L XL

MATERIAL

Styrene

TYPE

Perforated

12 per bag
Assorted bag: 16 trays



Solid

ITEM	EDENTULOUS
0921431	#1 Upper
0921432	#2 Upper
0921433	#3 Upper
0921434	#4 Upper
0921435	#5 Lower
0921436	#6 Lower
0921437	#7 Lower
0921438	#8 Lower
0921335	Set 104 contains 8 trays: 4U, 4L

ITEM	REGULAR
0921441	#1 Upper
0921442	#2 Upper
0921492	#3 Upper
0921493	#4 Upper
0921443	#5 Upper
0921444	#1 Lower
0921445	#2 Lower
0921494	#3 Lower
0921495	#4 Lower
0921446	#5 Lower
0921440	Set 105 contains 5 trays: 5U, 5L
0921186	#0 Upper (XL)
0921239	#0 Lower (XL)

ITEM	PEDO/ORTHO
0921720	#15 Upper
0921721	#16 Upper
0921722	#17 Upper
0921723	#18 Upper
0921724	#15 Lower
0921725	#16 Lower
0921726	#17 Lower
0921727	#18 Lower
0921719	Set 107 contains 8 trays: 4U, 4L



STAINLESS STEEL

STAINLESS STEEL TRAYS

100% Stainless Steel Trays

DESCRIPTION

Stainless Steel Impression Trays are created for use with alginate or any other type of impression material where detailed reproduction and dimensionally stable results are vital. They lock impression material in the tray to prevent distortion when removing it from the mouth. Uniform thickness gives superior rigidity for consistent results.

SIZES

S M L

MATERIAL

100% Stainless Steel

TYPE

Solid
Perforated

Edentulous: Sets of 8
Adult: Sets of 10



Perforated



#1 UPPER

SIZE
Small

ITEM
#0921600

(W) 2.67", 67.82 mm
(H) 2.02", 51.31 mm
(D) 0.61", 15.49 mm

#2 LOWER

SIZE
Small

ITEM
#0921599

(W) 2.26", 57.404 mm
(H) 1.92", 48.77 mm
(D) 0.48", 12.19 mm

#3 UPPER

SIZE
Medium

ITEM
#0921700

(W) 2.86", 72.64 mm
(H) 2.13", 54.10 mm
(D) 0.65", 16.51 mm

#4 LOWER

SIZE
Medium

ITEM
#0921701

(W) 2.74", 69.59 mm
(H) 2.14", 54.35 mm
(D) 0.56", 14.22 mm

#5 UPPER

SIZE
Large

ITEM
#0921702

(W) 3.01", 76.45 mm
(H) 2.32", 58.93 mm
(D) 0.71", 18.03 mm

#6 LOWER

SIZE
Large

ITEM
#0921703

(W) 2.86", 72.64 mm
(H) 2.15", 54.61 mm
(D) 0.56", 14.22 mm

#7 UPPER

SIZE
X-Large

ITEM
#0921706

(W) 3.07", 77.98 mm
(H) 2.62", 66.55 mm
(D) 0.71", 18.03 mm

#8 LOWER

SIZE
X-Large

ITEM
0921707

(W) 3.02", 76.71 mm
(H) 2.26", 57.40 mm
(D) 0.56", 14.22 mm

VARIETY BAG

SIZE
Contains 2 of
each size (16/bag)

ITEM
#0921704



BIG BITE TRAYS

Bite Registration Trays

DESCRIPTION

Big Bite Trays are economical, disposable, one piece plastic frames with a gauze insert to be used when recording centric occlusion. Big Bite Trays can be used on either side covering a quadrant or for a bi-lateral registration by clipping together right and left trays. The thin gauze mesh does not interfere with closure and ensures an accurate bite.

SIZES

One Size


MATERIALS

Plastic
Gauze Mesh

TYPE

Mesh

50 per box



Perforated Partial



7



8



9



10

#7 UR, LL

SIZE
Partial

ITEM
#0921370

(W) 2.24", 56.9 mm
(H) 2.29", 58.3 mm
(D) 0.79", 20.1 mm

#8 UL, LR

SIZE
Partial

ITEM
#0921371

(W) 2.15", 54.5 mm
(H) 2.33", 59.3 mm
(D) 0.86", 21.8 mm

#9 LOWER

SIZE
Anterior

ITEM
#0921372

(W) 2.56", 65.0 mm
(H) 1.50", 38.1 mm
(D) 0.75", 19.0 mm

#10 UPPER

SIZE
Anterior

ITEM
#0921373

(W) 2.52", 64.1 mm
(H) 1.47", 37.2 mm
(D) 0.76", 19.2 mm



50 PER BOX (25 L, 25 R)

SIZE
One Size

ITEM
#0921841

(W) 1.12", 28.5 mm
(H) 1.75", 44.54 mm
(D) 0.37", 9.5 mm



BOTRAYS®
Disposable Impression Trays

DESCRIPTION


BoTrays® Disposable Impression Trays are a sturdy, simple, and economical choice when consistent impressions are needed. These disposable trays are made from rigid, high-impact FDA approved styrene and organic FDA approved dye. Alginate "grippers" and slotted walls allow the BoTrays® to provide consistent, unsurpassed impressions with reduced gag effect.

SIZES
 S M L

MATERIAL
 Styrene

TYPE
 Perforated

12 per bag



Perforated



#1 UPPER

SIZE
 Large

ITEM
 #0921364

(W) 3.06", 77.6 mm
(H) 2.41", 61.2 mm
(D) 0.88", 22.4 mm

#2 LOWER

SIZE
 Large

ITEM
 #0921365

(W) 3.02", 76.7 mm
(H) 2.30", 58.3 mm
(D) 0.83", 21.1 mm

#3 UPPER

SIZE
 Medium

ITEM
 #0921366

(W) 2.91", 73.9 mm
(H) 2.21", 56.2 mm
(D) 0.83", 21.3 mm

#4 LOWER

SIZE
 Medium

ITEM
 #0921367

(W) 2.93", 74.3 mm
(H) 2.18", 55.2 mm
(D) 0.79", 20.0 mm

#5 UPPER

SIZE
 Small

ITEM
 #0921368

(W) 2.66", 67.9 mm
(H) 2.00", 50.9 mm
(D) 0.85", 21.9 mm

#6 LOWER

SIZE
 Small

ITEM
 #0921369

(W) 2.78", 70.6 mm
(H) 2.06", 52.3 mm
(D) 0.77", 19.6 mm