


EndoPilot®


Cordless Endodontic Handpiece

Based on dentists' operation habits, Endopilot provides cordless endodontic hadpiece in ergonomic design and small size. Therefore dentists can get better sight and more convenient operation.



Quick Start

Only one button can achieve root canal length measurement and preparation. The doctor can set the most commonly used parameters to this function, such as speed, torque, rotation mode.




Portable, Smart & Cordless

Cordless endodontic handpiece design. Safe and reliable medical equipment to let doctors use more handy.



Multi - position type head

It allows direction of contra angle to be changed easily. 6 angles available both hands are applicable. And it's convenient to use with push button.



Auto Reverse

The auto reverse function is activated depending on the Comparison between load and set torque. A wide variety of automatic functions are available.

Technical Data


Type:	Portable
Memory:	5 memory setting + quick start memory
Rotation:	Clockwise / Anti-clockwise / Reciprocating
Torque Control:	Auto reverse / Auto stop
Battery:	Rechargeable / Replaceable
Charging:	Charging base / USB
Display:	LED
Available Contra Angle:	10:1 16:1 20:1

Order Information

Part and its accessories:




Code: **AL2030**
Model: **EndoPilot**



Code: **ALP101**
Model: **16:1 contra angle**



Code: **ALP119**
Model: **Silicone Sleeve**



Code: **ALP111**
Model: **Spray nozzle**



Code: **ALP107**
Model: **Adaptor**



Code: **ALP104**
Model: **Battery**



Code: **ALP108**
Model: **Charging wire**



Code: **ALP121**
Model: **Charging base**

Beyes[®]

The Next Standard[®]

Beyes Dental Canada Inc.
30 Casebridge Court
Toronto, Ontario, M1B 3M5
Canada
Tel: 1-416-281-1888
Email: info@beyes.ca
Website: www.beyes.ca

Your Dealer for Beyes Products



Printed in Canada.
ENB002R1

© Beyes Canada Inc 2017. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.