



Disposable Plastic

## Multicolored Tips AIR/WATER SYRINGE TIPS

Multicolored Disposable Plastic Air Water Syringe Tips have a special locking grooves and are interchangeable with metal syringe tips, no adapter required. Quad air-core design with separate air and water channels to eliminate crossover. Fits all standard A-W syringes.

### FEATURES & BENEFITS

- Easy to use - no adapter needed, saves chair side time
- Superior air/water spray
- Unique quad aircore design ensures dry air when you need it
- Instant change from water to dry air
- Completely beveled - Smooth edge will not damage air water syringe O-ring
- Made from FDA approved plastics

The MARK3 Line  
was created to provide  
Premium Products  
at an Affordable Price

#1202 250 tips per bag. 60 bags per case.

#1203 1500 tips per bag. 8 bags per case.



Authorized for sale in Canada

*Compares to: Crystal Tips, Flashtips and Sparkle*