

SURGICAL CONTENTS

Aspirating Syringes	32
Anesthesia Instruments.....	33
Mouth Gags	34
Retractors	35
Scalpel Handles.....	35
Periosteal Elevators	36
Elevators	37-41
Cross-Bar Elevators	42-43
Root Tip Instruments	43-44
Fragment Forceps	44
Extracting Forceps	45-65
Rongeurs	66-69
Bone Curettes	70-71
Bone Files & Chisels.....	71-73
Mallets	74
Furnace Tongs	75
Dressing and Tissue Forceps	76-78
Hemostats	79-81
Needle Holders.....	82-92
Tungsten Carbide Needle Holders	88-92
Scissors.....	93-108
Tungsten Carbide Scissors	106-108
Towel Clamps	109



Aspirating Syringes

SURGICAL



CW-type Hub

Aspirating Syringe

07-802

CW-type Metric Hub

07-807



A-type Hub

Aspirating Syringe

07-801



CW-type

Aspirating Syringe

07-800

CW-Metric

Aspirating Syringe


07-806



A-type

Aspirating Syringe

07-805



Pool Suction Tube

8.75" Curved 23 French

07-880 Tips



Yankauer Suction Tube 11.25" *

07-885



Young Tongue Holder 6.25"

With Rubber Tips

24-4500

SURGICAL

* Not actual size



Silicon Replacement Tips
For Molt Mouth Gags

07-905



Molt Mouth Gag
Child

07-900



Molt Mouth Gag
Adult

07-910



#2 Austin Cheek Retractor

07-500



U/Minnesota Cheek Retractor

07-505



#3
For Blades 10, 11
12, 12B, 15, 15C
Scalpel Handle

19-3000



#4
For Blades 20-23
Scalpel Handle

19-4000



#7
Scalpel Handle

19-7000

Periosteal Elevators

SURGICAL



Freer
Periosteal Elevator

04-100



Molt #9
Periosteal Elevator

04-120



Goldman-Fox #14
Periosteal Elevator

04-130



OHL
Periosteal Elevator

04-140



Seldin #23
Periosteal Elevator

04-160



West #2
Periosteal Elevator

04-180

Our elevators feature all stainless steel construction for greater durability. Why is all stainless construction so important? Because it allows us to argon weld our tips to the handle – not simply press fit them as is done in brass handle elevators. This welding process gives us unsurpassed strength at the critical tip handle juncture, providing a stronger, more durable instrument. No chrome finish to wear or flake off, no press fit tip to work loose. Each perfectly forged tip undergoes a heat treatment in our high-tech, computerized vacuum furnace to assure a working end that will be forged to the optimum hardness without the worry of brittleness.



Apical #301
Elevator

04-401



Apical #302
Elevator

04-402



Apical #303
Elevator

04-403



Apical #304
Elevator

04-404



Bayonet #1
Elevator

04-201



Bayonet #2
Elevator

04-202

Elevators

SURGICAL



Crane Pick #41
Elevator

04-341



Cryer #45
Elevator Small Right

04-945



Cryer #26
Elevator Medium

04-926



Cryer #32
Elevator Large

04-332



Seldin #1L
Elevator

04-520



Cryer #44
Elevator Small Left

04-944



Cryer #25
Elevator Medium

04-925



Cryer #31
Elevator Large

04-331



Seldin #1R
Elevator

04-540



Seldin #3
Elevator
04-903



Seldin #4L
Elevator
04-560



Seldin #4R
Elevator
04-580



Heidbrink #1
Elevator
04-901



Heidbrink #2
Elevator
04-902



Heidbrink #40
Elevator
04-940



Heidbrink #41
Elevator
04-941



Heidbrink #80
Elevator
04-980

Elevators

SURGICAL



Schmeckebeer #81
Elevator

04-961



Howard #2
Elevator

04-982



Schmeckebeer #4
Elevator

04-964



Miller #71
Elevator

04-971



Miller #73
Elevator

04-973



Howard #1
Elevator

04-981



Howard #3
Elevator

04-983



Schmeckebeer #5
Elevator

04-965



Miller #72
Elevator

04-972



Miller #74
Elevator

04-974



Spear #36
Elevator

04-336



Universal #50
Elevator

04-350



General #190
Elevator

04-990



Gouge #34
Elevator

04-334



Gouge #46
Elevator

04-346



Stout #11A
Elevator

04-914



Universal #51
Elevator

04-351



General #191
Elevator

04-991



Gouge #34S
Elevator

04-335



Stout #11
Elevator


04-912

Cross-Bar Elevators

SURGICAL



#11L
Cross-Bar Elevator
04-953



#11R
Cross-Bar Elevator
04-954



#12L
Cross-Bar Elevator
04-955



#12R
Cross-Bar Elevator
04-956



#14L
Cross-Bar Elevator
04-957



#14R
Cross-Bar Elevator
04-958

Cross-Bar Elevators



Potts #1
Cross-Bar Elevator

04-591



Potts #2
Cross-Bar Elevator

04-592

Root Tip Picks & Elevators



Heidbrink #2/3
Root Tip Pick

04-200



Heidbrink #1
Root Tip Pick

04-210



Heidbrink #2
Root Tip Pick

04-220



Heidbrink #3
Root Tip Pick

04-230

SURGICAL

Root Tip Picks & Elevators

Fragment Forceps

SURGICAL



Root Tip #78

Elevator

04-378



Root Tip #79

Elevator

04-379



Root Tip #80

Elevator

04-380



Root Tip #341

Elevator

04-934



Stieglitz 90° Tip

Fragment Forceps 5"


04-389



Stieglitz Straight Tip

Fragment Forceps 5"

04-390




Stieglitz 45° Tip

Stieglitz Fragment Forceps 5"

04-385

J&J Instruments Extracting Forceps are the ideal complement to your surgical instrument inventory. Expertly crafted and carefully detailed, every pattern is correctly designed and every surface has been ground, polished, and/or beveled to perfection. The jaws have been forged to allow for anatomical nuances, facilitate grasping and maximize torque. Ten separate inspection stations assure that material and workmanship exceed industry standards. You can pay more for your extracting forceps, but why would you?

Most forceps can be serrated. Ask your dealer.



- 10 separate inspection stations
- 63 different patterns

#1 STD Standard
Upper Incisor Cuspid

05-001

Extracting Forceps

SURGICAL



90°
Silver Point Forceps 5"

04-395



45°
Fragment Forceps 5"

04-398



#3F
Woodward Lower 1st &
2nd Molar Universal Horn
Beak Straight Handle

05-030



#3FH
Woodward Lower 1st & 2nd Molar
Universal Horn Beak Hook Handle

05-034



#5
Lower 3rd Molar
Universal

05-050



#6
Splitting & Separating

05-060

Extracting Forceps

SURGICAL



#10H
Upper Molar
Universal Hook Handle

05-100



#10S
Upper Molar
Universal Straight Handle

05-105



#15
Lower 1st & 2nd Molar
Universal Hook Handle

05-150



#16
Lower 1st & 2nd Molar
Universal Hook Handle

05-160



#16S
Child Lower 1st & 2nd Molar
Universal Hook Handle

05-165



#17
Lower 1st & 2nd Molar
Universal Straight Handle

05-170

SURGICAL

Extracting Forceps

SURGICAL



#18L
Upper 1st & 2nd Molar Left
Hook Handle

05-180

#18R
Upper 1st & 2nd Molar Right
Hook Handle

05-185

#23
Lower 1st & 2nd Molar
Horn Beak Straight Handle

05-230



#24
Upper Molar Universal
Hook Handle

05-240

#32
Upper Universal
Bayonet-Bicuspid, Molar

05-320

#32A
Upper Universal
Bayonet-Bicuspid, Molar

05-321

SURGICAL

Extracting Forceps

SURGICAL



#53L

Upper 1st & 2nd Molar
Bayonet Left

05-530



#53R

Upper 1st & 2nd Molar
Bayonet Right

05-535



#62

Upper or Lower Universal
Incisor, Bicuspids, Deciduous

05-620



SURGICAL

#65
Upper Bicuspid Incisor
Root Bayonet

05-650

#69
Upper or Lower
Fragment or Small Root


05-690

#85
Lower Molars Cuspids
Bicuspids

05-850

Extracting Forceps

SURGICAL



#85A

Lower Molars Cuspids
Bicuspidis

05-851

#88L

Upper Left
1st & 2nd Molar

05-880

#88R

Upper Right
1st & 2nd Molar

05-885



#89
Upper Right Molar

05-890



#90
Upper Left Molar

05-900



#99A
Upper Incisor, Bicuspid,
Cuspid

05-991

SURGICAL

Extracting Forceps

SURGICAL



#99C

Upper Incisor, Bicuspid,
Cuspid

05-993



#101

Upper or Lower Universal
Bicuspid Deciduous

06-101



#103

Universal Bicuspid
Incisor Root

06-103



#150XS
Universal Bicuspid, Root Upper,
Pediatric with Spring

06-158



#151XS
Universal Incisor, Bicuspid Lower,
Pediatric with Spring

06-159

Extracting Forceps

SURGICAL



#150
Upper Incisor
Bicuspid Root Universal

06-150



#150 Serrated
Upper Incisor
Bicuspid Root Universal

06-156



#150S Child
Universal Bicuspid
Root Upper

06-154



#150A
Universal Incisor
Bicuspid Upper

06-152



SURGICAL

#151
Universal Bicuspid
Root Lower

06-151

#151 Serrated
Universal Bicuspid
Root Lower

06-157

#151A
Universal Incisor
Bicuspid Lower

06-153

#151S Child
Universal Incisor
Bicuspid Lower

06-155

Extracting Forceps

SURGICAL



#201
Extracting Forcep

06-201



#203
Universal, Lower Bicuspid,
Incisor Root

06-203



#210
Universal, Upper Molar

06-210



#210S
Universal, Upper Molar

06-211



#217
Universal
1st & 2nd Lower Molars

06-217



#222
Universal
Lower 3rd Molar

06-222

SURGICAL

Extracting Forceps

SURGICAL



#286

Upper Bicuspid
Incisor Root, Bayonet

06-286



#287

Universal
Lower 1st & 2nd Molar

06-287



#288

Lower, Bicuspid, Cuspid

06-288



English #73

Lower Molar

05-730



English #74

Lower Root

05-740



English #74N

Lower Root

05-745

SURGICAL

Extracting Forceps



English #13
Lower Bicuspid

05-130



English #22
Lower Molar

05-220



English #33
Lower Root

05-330

SURGICAL



#MD1

Universal Upper, Central,
Lateral, Cuspid, Bicuspid

06-321

#MD2

Universal Upper
1st, 2nd, 3rd Molar

06-322

#MD3

Universal Lower
Incisor, Cuspid, Bicuspid, Root

06-323

#MD4

Universal Lower
1st, 2nd, 3rd Molar


06-324

SURGICAL

Spring Kit (Not Shown)

Rongeur

11-010



Friedman 5.5"
30° Tip Rongeur

11-110



Mini Friedman 4.5"
30° Tip Rongeur

11-115



30°
11-030



45°
11-045



90°
11-090



Mini Blumenthal 4.5" Delicate

30° Tip Rongeur

11-020



Blumenthal 6"

30° Tip Rongeur

11-030


45° Tip Rongeur

11-045

90° Tip Rongeur

11-090

SURGICAL



#4
Rongeur 5.5"
11-140




#4A
Rongeur 6.5"
11-145



#5
Rongeur 6.5" Side Cutting

11-250



#10 (1A)
Rongeur 5.5"

11-310

Bone Currettes

SURGICAL



Lucas #75
Curette
02-275



Lucas #85
Curette
02-285



Lucas #86
Curette
02-286



Lucas #87
Curette
02-287



Lucas #88
Curette
02-288

Bone Curettes

Bone Files



Molt #1
Curette

02-291



Molt #2
Curette

02-292



Molt #4
Curette

02-294



Molt #5L
Curette

02-295



Molt #6R
Curette

02-296



#12A
Bone File Cross-Serrated

02-100

SURGICAL

Bone Files

SURGICAL



#10
Bone File

02-110




Howard #12
Bone File

02-120



Miller #21
Bone File

02-210



Miller #64
Bone File

02-364



Miller #45
Bone File

02-450



Coupland #1
Bone Chisel

02-010



Coupland #2
Bone Chisel

02-020



Gardner #51
Bone Chisel

02-310



Gardner #52
Bone Chisel

02-320



Wakefield #1W
Bone Chisel

02-370



Wakefield #2W
Bone Chisel

02-375

SURGICAL

Mallets

SURGICAL



#49
Mallet 7.5" Lead-Faced

02-350



#49
Mallet 7.5" Steel-Faced

02-352



Mead Mallet
7.5" Rubber Faced

02-360



Replacement Facing
Mead Mallet Replacement

02-365



10" *
Furnace Tongs Serrated

35-0100



12" *
Furnace Tongs Serrated

35-0120

SURGICAL

* Not actual size

Dressing and Tissue Forceps

SURGICAL



5"
Dressing FCP Serrated

35-0050



5.5"
Dressing FCP Serrated

35-0055



6"
Dressing FCP Serrated

35-0060



8"
Dressing FCP Serrated

35-0080



5"
Tissue FCP 1x2

35-1250



5.5"
Tissue FCP 1x2

35-1255



6"
Tissue FCP 1x2

35-1260



Adson 4.75"
Dressing FCP Serrated

35-5010

Serrated
Tissue FCP
35-5010


SURGICAL

Dressing and Tissue Forceps

SURGICAL



1 x 2
Tissue FCP
35-5020



24-3010 **24-3050** **24-3090**

7 x 7
Tissue FCP
35-5090




Adson 4.75" 1x2
Tissue FCP
35-5020



Adson-Brown 4.75" 7x7
Tissue FCP
35-5090



Allis 6" 5x6
Tissue FCP
24-3050



Allis 7.5" 5x6
Tissue FCP
24-3090


Allis 6" 4x5
Tissue FCP
24-3010

Examine the JJI hemostats and needle holders carefully. You'll notice that all the box lock joints are meticulously beveled. This not only gives each instrument a finished quality look, but also reduces the possibility of snagging suture or exam gloves. Notice too, as the jaws begin to touch so do the ratchets. This results in all JJI forceps and needle holders having a consistent balance, feel, and force-to-closure ratio.



Mosquito 3.5"
Curved FCP
24-1015

Mosquito 3.5"
Straight FCP
24-1010



Mosquito 5"
Curved FCP
24-1055

Mosquito 5"
Straight FCP
24-1050

Hemostats

Rochester Pean 5.5"

Straight FCP

24-1300

Rochester Pean 5.5"

Curved FCP

24-1305

Rochester Pean 6.25"

Curved FCP

24-1325

Rochester Pean 6.25"

Straight FCP

24-1320



Kelly 5.5"

Curved FCP

24-1135



Kelly 5.5"

Straight FCP

24-1130

Rankin-Kelly 6.25"

Curved FCP

24-1175



Rankin-Kelly 6.25"

Straight FCP

24-1170



Crile 5.5"
Curved FCP
24-1245

Crile 5.5"
Straight FCP
24-1240



Rankin-Crile 6.25"
Curved FCP
24-1285


Rankin-Crile 6.25"
Straight FCP
24-1280

SURGICAL



Derf 4.75"
Needle Holder Serrated

25-6000



Webster 5"
Needle Holder Smooth

25-6050



Halsey 5"
Needle Holder Smooth


25-6070



Baumgartner 5.25"
Needle Holder

25-6210

Needle Holders



SURGICAL



Mayo-Hegar 5"

Needle Holder

25-6290

Mayo-Hegar 5.5"

Needle Holder

25-6300

Mayo-Hegar 6"

Needle Holder

25-6360

Mayo-Hegar 7"

Needle Holder

25-6370



Crile-Wood 6"

Needle Holder

25-6420



Collier 5"
Fenestrated Jaws
25-7100


Olsen-Hegar 5.5"
Needle Holder
25-6530

Olsen-Hegar 6.5"
Needle Holder
25-6540

Abbey 5"
Needle Holder
25-7180

Needle Holders

SURGICAL



Boynton 5"
Needle Holder

25-7550



Mathieu 4.75" Delicate
Needle Holder


25-7680



25-7650



25-7655



Mathieu 5.5" Regular
Needle Holder

25-7650

Mathieu 5.5" Delicate
Needle Holder

25-7655



Castroviejo 5.75" No Lock
Needle Holder

25-7800



Castroviejo 5.75" W/ Lock
Needle Holder

25-7850

Tungsten Carbide Inserted Needle Holders

Gator-Grip™ instruments by J&J Instruments are our top of the line needle holders. From their sleek, gold-plated handles to the German made tungsten-carbide tipped inserts, these needle holders have the feel, quality and durability to please the most demanding surgeon. The tungsten inserts are affixed by a copper brazing process and heat treated in our computerized vacuum furnace. This is done to assure the strongest possible bond between the insert and body of the needle holder. This bond will remain intact for a minimum of 10 years – we guarantee it.

SURGICAL




Mayo-Hegar 5.5" Gator Grip™

Tungsten Carbide Needle Holder

65-6300

Tungsten Carbide Inserted Needle Holders




Olsen-Hegar 6.25"
Tungsten Carbide Needle Holder
65-6540


Olsen-Hegar 7.5"
Tungsten Carbide Needle Holder
65-6550



Mayo-Hegar 7" Gator Grip™
Tungsten Carbide Needle Holder
65-6370



Derf 4.75" Gator Grip™
Tungsten Carbide Needle Holder
Serrated Jaws
65-6000



Olsen-Hegar 5.5"
Tungsten Carbide Needle Holder
65-6530



Mayo-Hegar 8" Gator Grip™
Tungsten Carbide Needle Holder
65-6380



Mayo-Hegar 6" Gator Grip™
Tungsten Carbide Needle Holder
65-6360

SURGICAL

Tungsten Carbide Inserted Needle Holders

SURGICAL



Crile-Wood 6" Gator Grip™
Tungsten Carbide Needle Holder

65-6420



Boynton 5" Gator Grip™
Tungsten Carbide Needle Holder

65-7550



Mathieu 5.5" Gator Grip™
Tungsten Carbide Needle Holder

65-7650



Mathieu 5.5" Gator Grip™
Delicate Carbide Needle Holder

65-7655



Mathieu 4.75" Gator Grip™
Delicate Carbide Needle Holder

65-7680

Tungsten Carbide Inserted Needle Holders


SURGICAL



Castroviejo 5.75" Gator Grip™

Tungsten Carbide Needle Holder

65-7850



Castroviejo 7" Gator Grip™

Tungsten Carbide Needle Holder

65-7875

There simply is no better scissor value in the dental industry today. Every J&J Instruments scissor is forged from German series 420B stainless steel. Each pattern has a certain delicate balance and feel, comfortable rings, easy cutting action and durability. If you come across a scissor to match our features ... compare the price. You'll understand.



Goldman-Fox 5"
Curved, 1 Serrated Blade Scissor


12-125



Dean 7"
Angled, 1 Serrated Blade Scissor

12-100

SURGICAL



Crown 4.5"
Straight Scissor

12-010



Crown 4.5"
Straight Serrated Scissor

12-030



Crown 4.5"
Curved Scissor

12-020



Crown 4.5"
Curved Serrated Scissor

12-040





Kelly 6.25"
Straight Scissor

12-215

Kelly 6.25"
Straight Scissor w/ 1 Serrated Blade

12-275



Kelly 6.25"
Angled Scissor

12-195

Kelly 6.25"
Angled Scissor w/ 1 Serrated Blade

12-255



Kelly 6.25"
Curved Scissor

12-235

Kelly 6.25"
Curved Scissor w/ 1 Serrated Blade

12-295

SURGICAL



Wagner 4.75"

Straight Scissor

12-318

Wagner 4.75"

Straight Scissor w/ 1 Serrated Blade

12-339



Wagner 4.75"

Angled Scissor

12-306

Wagner 4.75"

Angled Scissor w/ 1 Serrated Blade

12-330



Wagner 4.75"

Curved Scissor

12-324

Wagner 4.75"

Curved Scissor w/ 1 Serrated Blade

12-345




Quimby 5"
Straight Scissor

12-420



Quimby 5"
Curved Scissor

12-440



5.5"
Straight Operating Scissor, Sharp/Sharp

22-1100


5.5"
Straight Operating Scissor, Sharp/Blunt

22-1110



Iris 4.5"
Straight Scissor

22-1500



Iris 4.5"
Curved Scissor

22-1510



La Grange 4.5"
Curved Scissor

22-1540



Castroviejo 4.75"
Curved Scissor

22-1595



Mayo 5.5"
Straight Scissor

22-1600



Mayo 5.5"
Curved Scissor

22-1605

Scissors (Suture)

SURGICAL



Spencer 3.5"
Stitch Scissor

23-3010



Littauer 4.5"
Stitch Scissor

23-3040



Shortbent 3.5"
Stitch Scissor

23-3030



Littauer 5.5"
Stitch Scissor

23-3060



6.25"
Extra Long Stitch Scissor

23-3070

Scissors (Suture)

Scissors (Wire Cutting)

SURGICAL



Northbent 4.75"
Curved Stitch Scissor

23-3085



4.75"
Angled Wire Cutting Scissor

23-4050

SURGICAL



Lister 5.5"
Bandage Scissor

23-4250



Lister 3.5"
Bandage Scissor

23-4200



Lister 4.5"
Bandage Scissor

23-4230



Lister 7.25"
Bandage Scissor

23-4270

Scissors (Bandage & Utility)

SURGICAL



Knowles 5.5"
Bandage Scissor

23-6200



7.5"
Utility Scissor Black

23-7070

7.5"
Utility Scissor Blue

23-7071

7.5"
Utility Scissor Gray

23-7072

7.5"
Utility Scissor Green

23-7073



8"
US Army Gauze Scissor

23-7570



Castroviejo 3.5"
Straight Scissor

22-1580



Castroviejo 3.5"
Angled Scissor

22-1585



Castroviejo 4.75"
Straight Scissor

22-1590

Scissors (Gator-Edge w/ Carbide Inserted Blades)

The Gator-Edge scissors are manufactured for the practitioner who demands the best. We start with the finest materials. From the German 420B stainless steel body of the scissor to the German TC inserts. Then every tungsten carbide inserted cutting edge is precision-honed to a durable, long-lasting and sharp surface. The edge complements the delicate feel and balance of the scissor to provide years of smooth and dependable use. Look for the gold rings. Feel the difference.

SURGICAL



Dean 7"
Angled, Carbide Scissor

62-100



Goldman-Fox 5"
Curved, Carbide Scissor

62-125

Scissors (Gator-Edge w/ Carbide Inserted Blades)

SURGICAL



Iris 4.5"
Straight, Carbide Scissor

62-1500




Iris 4.5"
Curved, Carbide Scissor

62-1510

Scissors (Gator-Edge w/ Carbide Inserted Blades)

SURGICAL



Mayo 5.5"
Straight, Carbide Scissor

62-1600



Mayo 5.5"
Curved, Carbide Scissor

62-1605



Kelly 6.25"
Angled, Carbide Scissor

62-195



Kelly 6.25"
Straight, Carbide Scissor

62-215



Kelly 6.25"
Curved, Carbide Scissor

62-235



Jones 3.5"
Towel Clamp
31-1190



Backhaus 3.5"
Towel Clamp
31-1000



Backhaus 5.25"
Towel Clamp
31-1010