


# Sonic 100

Aurelia® Sonic® Nitrile Indigo Blue  
Nitrile Examination Gloves



Indigo Blue Color  
NEW 100/box pack

- New ultra comfort formulation
- 2.2 mil thickness (palm measurement) for improved tactile sensitivity
- Finger textured finish
- 100 gloves per box
- Soothing indigo blue color



**AURELIA**  
Sign of Comfort™

[www.aureliagloves.com](http://www.aureliagloves.com) | 1-877-AURELIA (1-877-287-3542)

# Sonic 100

Our revolutionary and popular ultra-thin formulation Sonic gloves now come in machine packed 100 count packs for more sonic-paced convenience and ease of ordering. Sonic 100 features the same renowned high tactility and fashionable lightweight comfort of Sonic but in a smaller count packaging to better complement today's more dynamic PPE environment. For more information please visit [www.aureliagloves.com](http://www.aureliagloves.com).



## Aurelia® Sonic 100™ Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	680
Small	31.1	680
Medium	31.0	680
Large	31.2	680
Extra Large	31.3	680
ASTM D6319	14 MPa minimum	500% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	620
Small	30.0	620
Medium	30.0	620
Large	30.0	620
Extra Large	30.0	620
ASTM D6319	14 MPa minimum	620% minimum

### DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>Extra Small</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Small</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Medium</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Large</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Extra Large</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum

### Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

### ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large*
9377B5	9377B6	9377B7	9377B8	9377B9

Authorized Aurelia Distributor:



**AURELIA**  
Sign of Comfort™

AURELIA Div. Supermax Inc.

1899 Sequoia Drive, Aurora, IL 60506 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855  
Web: [www.aureliagloves.com](http://www.aureliagloves.com) | Email: [info@aureliagloves.com](mailto:info@aureliagloves.com) | [facebook.com/aureliagloves](https://facebook.com/aureliagloves)