



## A-Dec Style

# Metal 3-Way Tips AIR/WATER SYRINGE TIPS

These A-dec style, autoclavable, air and water syringe tips have three locking grooves that hold the tip in place when properly installed and prevent ejection of the tip during use. Stainless steel syringe tips are specifically designed for A-dec type syringes.

### FEATURES & BENEFITS

- Locking grooves hold tips in place
- Stainless steel specifically designed for A-Dec type syringes
- Eliminates possibility of cross contamination

The MARK3 Line  
was created to provide  
Premium Products  
at an Affordable Price

**#1401** 5 tips per bag

*Compares to: Safe-Tips EZ*



Authorized for sale in Canada