

MARK3™
THE SMART ALTERNATIVE



Vinyl Polysiloxane • Super Hydrophilic VPS 380 Impression Material

FEATURES & BENEFITS

- New delivery system provides an economic way to take accurate impressions. For use in any dynamic mixing machine. Creates flawless, void free impressions in seconds.
- Excellent dimensional stability decreases distortions
- Regular Set working times
- No odor or irritation allowing greater patient comfort
- Pleasant mint flavor
- Made in the USA



REFILL KIT - 2 x 380ml cartridges per box

#3025Heavy Body - Regular Set

#3026Monophase - Regular Set

Work Time: 2.30 minutes


Intra Oral Set Time: 1 minute

Compare to Genie VPS Magic Mix

Accessories

Dynamic Yellow Mixing Tips

order# 1419
50 per bag



Bayonet Replacement Rings

order# 1420
2 per bag

