



PRE-BENT

Applicator Tips

MARK3 Pre-Bent Applicator Tips are blunt end, pre-bent Luer-Lock dispensing needle tips for fast, simple placement of medicaments, etchants, resins or flowable composites.

FEATURES & BENEFITS

- Stainless steel construction
- Angled cannula epoxy sealed to a double helix luer lock polypropylene hub
- Polished burr-free finish for smooth fluid dispensing
- Disposable blunt end
- For use with etchants, resins, and flowable composites

The MARK3 Line
was created to provide
High Quality Products
at an Affordable Price

- #2270** Pink Tips - 18 gauge
- #2280** Yellow Tips - 20 gauge
- #2290** Black Tips - 19 gauge
- #2300** Grey Tips - 22 gauge
- #2310** Blue Tips - 25 gauge

100 per pack



AUTHORIZED FOR SALE IN CANADA