



**GLUMA® Desensitizer & GLUMA® Desensitizer PowerGel**  
Under every restoration.

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

# GLUMA® Desensitizer

## GLUMA® Desensitizer PowerGel

Two solutions to cover all of your sensitivity needs



### GLUMA Desensitizer

GLUMA Desensitizer has been in the market for more than 20 years and has been used in more than 50 million restorations worldwide to stop and prevent hypersensitivity. One drop of Gluma is all you need and the result is fast and effective. No mixing, no curing and no repetitive steps.




### GLUMA Desensitizer PowerGel

With the same proven efficacy as GLUMA Desensitizer, GLUMA Desensitizer PowerGel is a unique one-step gel formula desensitizer that allows for greater control and accuracy during application. In addition, its unique green color indicator provides visual ease of placement and rinse-clean capability.

### How it works:

- **Desensitize:** GLUMA is the only desensitizer proven to penetrate exposed dentinal tubules up to 200µ. The result is the formation of multiple layers of protein septas that prevent inter-tubular fluid movements due to osmotic changes.
- **Inhibit bacteria:** GLUMA provides a hermetic seal that acts as a microbial barrier, inhibiting bacterial growth.
- **Re-wet:** GLUMA resurrects collapsed collagenous fibers, improving the bond strength of many adhesives.



Dentinal tubule before application of GLUMA



Protein septas following application of GLUMA

| GLUMA® Desensitizer                         | Order Number |
|---|--------------|
| GLUMA Desensitizer 5ml                      | 65872354     |
| GLUMA Desensitizer Clinic Pack 3 x 5ml      | 66018221     |
| GLUMA Desensitizer Single Dose 40 x 0.075ml | 66001854     |
| Application Brushes 2 Pack; 50 Pkg          | 50025123     |
| Brush Handle 1 Pack                         | 50025125     |

| GLUMA® Desensitizer PowerGel       | Order Number |
|------------------------------------|--------------|
| GLUMA PowerGel Syringe Kit 4 x 1gm | 66043451     |
| GLUMA PowerGel Brush Cannulas 60   | 66044334     |

FIRST  
DESENSITIZER  
ON THE MARKET  
1987



DENTAL ADVISOR  
1997 • 2000  
2008 (Top Desensitizer)



REALITY'S CHOICE  
FIVE STAR AWARD  
2001 • 2002 • 2003  
2005 • 2006 • 2007  
2009 • 2011



TOWNIE CHOICE  
AWARD  
2003 • 2004 • 2005  
2006 • 2007 • 2008  
2009 • 2010 • 2011

CLINICIAN'S  
REPORT  
2009

Top Effective Product for Class II

GLUMA POWER GEL



DENTAL ADVISOR  
2012



CPS FIVE  
DIAMOND AWARD  
2011



REALITY'S CHOICE  
FIVE STAR AWARD  
2013