

# CANAVIEW<sup>®</sup>

## Intraoral Camera


An excellent communication tool that improves understanding and interaction with patients. It provides visual evidence that you need to educate patients and make more accurate diagnosis. Boasting a true autofocus and intuitive camera lighting system, the dental digital photography camera requires no manual adjustments and is easily shared between different examination rooms and chairs.

- High Resolution images with 5MP
- Ultra-wide and non-spherical lens
- Innovative Anti-fog design
- Direct USB connection to computer
- Super light weight body
- Responsive capture button
- Complimentary Image management software
- Compatible with major imaging software:  
*AbleDent, Aptyx, Adstra, CaptureLink, CliniView, Cleardent, DentiMax, Digora, Dexis, Eaglesoft, Owandy, ProfSuni, and more...*



## Key Specifications

Focus Type	Auto
Focus Range	4mm - 60mm
Connection Type	USB to computer
Image Quality	5 million pixels
Resolution	1600 x 1200
Camera Length	20cm
Weight	35g
Cable Length	2m
Power	DC 5V+/-0.5V(1.5A)
Operating System	Windows 10/8/7/VXP/Vista (64/32 bit)
Warranty	12 months



## Order Information

Code: EV2060    Model: Canaview



**Beyes**<sup>®</sup>  
Partnership for Success

Beyes Dental Canada Inc.  
30 Casebridge Court  
Toronto, Ontario, M1B 3M5  
Canada  
Tel: 1-416-281-1888  
Email: info@beyes.ca  
Website: www.beyes.ca

Your Dealer for Beyes Products

Scan the QR code for  
more information on CanaView



Printed in Canada  
File #: ENB007R2

© Beyes Canada Inc. 2018. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance  
with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under Canadian  
Medical Devices Conformity Assessment System (CMDCAS).