

MARK3™
THE SMART ALTERNATIVE



10:1 Automix

Temporary Crown & Bridge

FEATURES & BENEFITS

- Strong and durable for short and long lasting provisionals
- Excellent mechanical and esthetic properties
- Quick and easy to use
- Easily removable
- Affordable price

Technical Data

- Dosage and application time: Approx. 45 seconds
- Setting time in mouth: Approx. 2-3 minutes (hard elastic phase, from beginning of mixing)
- Setting time: Approx. 6-7 minutes (until setting time is completed)
- Mixing ratio: 10:1

4101.....Shade A1

4102.....Shade A2

4103.....Shade A3

1 Cartridge of 50ml plus 10 Mixing Tips

Used with...

Dispensing
Gun - 10:1

order# 110411

Type 50

*Suggested Retail
Price: \$49.99*



**Compare to Luxatemp, Protemp,
Integrity and Temphase**